

## ABS\_Bio™ Wound Healing Assay Kit (Cat# K040-24; 24 assays; store kit at 4°C)

### Introduction

The wound healing assay has been widely adapted and modified to study the effects of a variety of experimental conditions, for instance, gene knockdown or chemical exposure, on mammalian cell migration and proliferation.

The ABS\_Bio™ Wound Healing Assay Kit provides a simple, sensitive, one-step colorimetric and fluorimetric assay for cell migration and proliferation. In this assay, the inserts create a wound field with a defined gap of 0.5mm and the “healing” of this gap by cell migration and growth towards the center of the gap is monitored and quantitated. Factors that alter the motility and/or growth of the cells can lead to increased or decreased rate of “healing” of the gap. The cell proliferation and migration rates can be determined using manual fixing and microscopic imaging. The kit is supplied with sufficient reagents for 24 tests in 24-well plate assay. It could easily be modified for a high-throughput assay.


### Kit Components (24 tests)

Fixation Solution: 40 mL    Light Stain Solution: 40 mL    Fluorescence Stain Solution(100x): 0.4 mL  
24-Well plate with 24 ready to use Culture-Inserts: 1

**Storage and Handling:** Store Fluorescence Stain Solution at -20°C, others components at 4°C. Shelf Life: 12 months after receipt. Warm up Reagents to room temperature before use.

### Procedure

1. Grow cells in DMEM supplemented with 10% FBS. Prepare cell suspension containing  $3-7 \times 10^5$  cells/mL in media.



2. Using either 24-Well plate with 24 ready to use Culture-Inserts (#S80241) or self-insertion Culture-Insert (#S80209) to process the wound healing assay.

3. If use self-insertion Culture-Insert (#S80209), using sterile tweezers, place the desired number of Culture-Inserts in the plate wells. **Note: Culture-Insert with special sticky and biocompatible surface at the bottom side works like a glue and avoids leaking.** It fit standard 12-Well plates, 6-Well plates or petri dishes. If use 24-Well plate with 24 ready to use Culture-Inserts (#S80241), directly start from step 4.

4. Apply 70  $\mu$ L cells into each insert well, they should reach ~80-90% confluence within overnight culture.

5. After appropriate cell attachment gently remove the Culture-Insert by using sterile tweezers. Grab a corner of the Culture-Insert.

6. If necessary, a washing step can help removing non-adherent cells or cell debris.




7. Replenish the well with fresh medium to start wound healing process.

**Note: Medium may contain ingredients of interest that you want to test, e.g., chemicals that inhibit/promote cell motility and/or proliferation.**

8. Grow cells for additional 24-48h (or the time required if different cells are used).

Monitor the wound closure with a light microscope or image assay (see below cell stain).

9. Measure the percent closure or the migration rate of the cells into the wound field.



## Calculation of Results


### Percent Closure:

1. Determine the surface area of the defined wound area. Total Surface Area = wide x length (mm)
2. Determine the surface area of the migrated cells in to the wound area. Migrated Cell Surface Area = wide of cell migration (mm) x 2 x length (mm)
3. Percent Closure (%) = Migrated Cell Surface Area / Total Surface Area x 100

### Migration Rate:

Determine the migration rate of cells into the defined wound area:

Migration Rate = length of cell migration (nm) / migration time (hr).



Cells light stain assay

### Cell stain Assay

1. Wash the cells twice with 1x PBS.
2. Fix the cells with 400-1000  $\mu$ L fixation solution for 15 min at RT
3. Aspirate and discard the solution. Wash the cells twice with 1x PBS.

**For light stain:** Stain the fixed cells with 400-1000  $\mu$ L Light Stain Solution for 15 min at RT, wash the cells twice with 1x PBS, discard washes and allow cells to dry at room temperature. Shoot pictures with light microscope.

**For fluorescence stain:** 100x dilute Fluorescence Stain Solution with PBS, add 400-1000  $\mu$ L diluted stain solution into each well for 15 min at RT, wash the cells twice with 1x PBS. Shoot pictures with Fluorescence microscope with 350nm/470nm filter.

The gap distance can be quantitatively evaluated using software such as Photoshop or ImageJ (<http://rsb.info.nih.gov/ij/download.html>)

## References

Liang, CC. et al. 2007, Nature Protocol 2:329-333  
Yarrow, JC. et al. 2004, BMC Biotechnology 4:21

## Related Products:

**24-Well plate with 24 ready to use Culture-Inserts for wound healing assay (#S80241)**  
**Culture-Insert for wound healing assay (#S80209)**

Trypan Blue dye (#C8039)

Cell Viability Detection Kit(WST-1 based; #K010-500)

Janus Green Cell Stain Solution (#C8070)

Wound Healing Assay Kit (#K040-25)

Cell Viability Detection Kit (LDH based; #K025-500)

Scratch Wound Healing Assay Kit (#K040-100)