

# **Human IgG Antigen ELISA Kit**

Catalog # HUIGGKT

Strip well format. Reagents for up to 96 tests. Rev: August 2017

#### INTENDED USE

This human immunoglobulin G antigen assay is intended for the quantitative determination of total human IgG antigen in serum, plasma, hybridoma cell supernatants, ascites or other biological fluids. The assay does not distinguish IgG subclasses. For research use only.

#### **BACKGROUND**

IgG is the most abundant immunoglobulin in serum and is predominately involved in the secondary immune response. The IgG subclasses are designated 1, 2, 3 and 4 based on their relative prevalence in human serum.

### **ASSAY PRINCIPLE**

Human IgG will bind to the affinity purified capture antibody coated on the microtiter plate. After appropriate washing steps, horseradish peroxidase labeled polyclonal anti-human IgG antibody binds to the captured protein. Excess antibody is washed away and TMB substrate is used for color development at 450nm. A standard calibration curve is prepared along with the samples to be measured using dilutions of human IgG. Color development is proportional to the concentration of IgG in the samples.

## STANDARD CALIBRATION

Human IgG standard provided is calibrated against the WHO International Standard for Immunoglobulins G, A and M, Human Serum distributed by NIBSC (67/086), South Mimms, Potters Bar, Hertfordshire, UK.

Lot 615L: 500 ng = 0.00525 IU

# **REAGENTS PROVIDED**

- •96-well antibody coated microtiter strip plate (removable wells 8x12) containing anti-human IgG antibody, blocked and dried.
- •10X Wash buffer: 1 bottle of 50ml
- Human IgG standard: 1 vial lyophilized standard
- Anti-human IgG primary antibody: 1 vial lyophilized polyclonal antibody
- Horseradish peroxidase-conjugated streptavidin: 1 vial concentrated HRP labeled streptavidin
- •TMB substrate solution: 1 bottle of 10ml solution

#### STORAGE AND STABILITY

Store all kit components at 4°C upon arrival. Return any unused microplate strips to the plate pouch with desiccant. Reconstituted standards and primary may be stored at -80°C for later use. Do not freeze-thaw the standard and primary antibody more than once. Store all other unused kit components at 4°C. This kit should not be used beyond the expiration date.

# OTHER REAGENTS AND SUPPLIES REQUIRED

- Microtiter plate shaker capable of 300 rpm uniform horizontally circular movement
- Manifold dispenser/aspirator or automated microplate washer
- Microplate reader capable of measuring absorbance at 450 nm
- Pipettes and Pipette tips
- Deionized or distilled water
- Polypropylene tubes for dilution of standard
- Paper towels or laboratory wipes
- •1N H<sub>2</sub>SO<sub>4</sub> or 1N HCl
- •Bovine Serum Albumin Fraction V (BSA)
- Tris(hydroxymethyl)aminomethane (Tris)
- Sodium Chloride (NaCl)

# **PRECAUTIONS**

- •FOR LABORATORY RESEARCH USE ONLY. NOT FOR DIAGNOSTIC USE.
- Do not mix any reagents or components of this kit with any reagents or components of any other kit. This kit is designed to work properly as provided.
- Always pour peroxidase substrate out of the bottle into a clean test tube. Do not pipette out of the bottle as contamination could result.
- Keep plate covered except when adding reagents, washing, or reading.
- DO NOT pipette reagents by mouth and avoid contact of reagents and specimens with skin.
- •DO NOT smoke, drink, or eat in areas where specimens or reagents are being handled.

### PREPARATION OF REAGENTS

- •TBS buffer: 0.1M Tris, 0.15M NaCl, pH 7.4 •Blocking buffer (BB): 3% BSA (w/v) in TBS
- •1X Wash buffer: Dilute 50ml of 10X wash buffer

concentrate with 450ml of deionized water

### **SAMPLE COLLECTION**

Collect plasma using EDTA or citrate as an anticoagulant. Centrifuge for 15 minutes at 1000xg within 30 minutes of collection. Assay immediately or aliquot and store at ≤ -20°C. Avoid repeated freeze-thaw cycles.

## **ASSAY PROCEDURE**

Perform assay at room temperature. Vigorously shake plate (300rpm) at each step of the assay.

## **Preparation of Standard**

Reconstitute standard by adding 1ml of blocking buffer directly to the vial and agitate gently to completely dissolve contents. This will result in a 500ng/ml standard solution.

Dilution table for preparation of human IgG standard:

IgG concentration (ng/ml)	Dilutions
200	600μl (BB) + 400μl (from vial)
100	500μl (BB) + 500μl (200ng/ml)
50	500µl (BB) + 500µl (100ng/ml)
20	600µl (BB) + 400µl (50ng/ml)
10	500µl (BB) + 500µl (20ng/ml)
5	500µl (BB) + 500µl (10ng/ml)
2	600µl (BB) + 400µl (5ng/ml)
1	500µl (BB) + 500µl (2ng/ml)
0.5	500µl (BB) + 500µl (1ng/ml)
0.2	600µl (BB) + 400µl (0.5ng/ml)
0	500µl (BB) Zero point to determine background

NOTE: DILUTIONS FOR THE STANDARD CURVE AND ZERO STANDARD MUST BE MADE AND APPLIED TO THE PLATE IMMEDIATELY.

# **Standard and Unknown Addition**

Remove microtiter plate from bag and add 100µl IgG standards (in duplicate) and unknowns to wells. Carefully record position of standards and unknowns. Shake plate at 300rpm for 30 minutes. Wash wells three times with 300µl wash buffer. Remove excess wash by gently tapping plate on paper towel or kimwipe.

NOTE: The assay measures total human IgG antigen in the 0.2-200 ng/ml range. If the unknown is thought to have IgG levels greater than 200ng/ml, dilutions may be made in blocking buffer before use. A 1:1,000,000-1:5,000,000 dilution for normal human plasma is suggested for best results.

## **Primary Antibody Addition**

Reconstitute primary antibody by adding 10ml of blocking buffer directly to the vial and agitate gently to completely dissolve contents. Add 100 $\mu$ l of this dilution to all wells. Shake plate at 300rpm for 30 minutes. Wash wells three times with 300 $\mu$ l wash buffer. Remove excess wash by gently tapping plate on paper towel or kimwipe.

# Streptavidin-HRP Addition

Briefly centrifuge vial before opening. Dilute 2.5 $\mu$ l of HRP conjugated streptavidin into 2.5 $\mu$ l blocking buffer to generate a 1:1,000 dilution. Add 0.1 $\mu$ l of 1:1,000 dilution to 9.9 $\mu$ l of blocking buffer to generate a 1:100,000 dilution. Add 100 $\mu$ l of the 1:100,000 dilution to all wells. Shake plate at 300 $\mu$ l wash buffer. Remove excess wash by gently tapping plate on paper towel or kimwipe.

#### **Substrate Incubation**

Add 100µl TMB substrate to all wells and shake plate for 2-10 minutes. Substrate will change from colorless to different strengths of blue. Quench reaction by adding 50µl of 1N  $H_2SO_4$  or HCl stop solution to all wells when samples are visually in the same range as the standards. Add stop solution to wells in the same order as substrate upon which color will change from blue to yellow. Mix thoroughly by gently shaking the plate.


# Measurement


Set the absorbance at 450nm in a microtiter plate spectrophotometer. Measure the absorbance in all wells at 450nm. Subtract zero point from all standards and unknowns to determine corrected absorbance ( $A_{450}$ ).

# **Calculation of Results**

Plot A<sub>450</sub> against the amount of IgG in the standards. Fit a straight line through the linear points of the standard curve using a linear fit procedure if unknowns appear on the linear portion of the standard curve. Alternatively, create a standard curve by analyzing the data using a software program capable of generating a four parameter logistic (4PL) curve fit. The amount of IgG in the unknowns can be determined from this curve. If samples have been diluted, the calculated concentration must be multiplied by the dilution factor.

A typical standard curve (EXAMPLE ONLY):





# **EXPECTED VALUES**

The concentration of IgG in normal human serum ranges from 5 to 12 mg/ml.

### PERFORMANCE CHARACTERISTICS

**Sensitivity:** The minimum detectable dose (MDD) was determined by adding two standard deviations to the mean optical density value of twenty zero standard replicates (range  $OD_{450}$ : 0.066-0.081) and calculating the corresponding concentration. The MDD was 0.071 ng/mL.

**Intra-assay Precision:** Three samples of known concentration were tested twenty times on one plate to assess intra-assay precision.

Sample		1	2	3	
n		20	20	20	
Mean (ng/ml)		0.729	2.31	49.5	
Standard Dev	iation	0.032	0.072	2.97	
CV (%)		4.36	3.13	6.00	

**Inter-assay Precision:** Three samples of known concentration were tested in ten independent assays to assess inter-assay precision.

1	Sample	1	2	3
	n	10	10	10
	Mean (ng/ml)	0.967	2.54	27.6
	Standard Deviation	0.050	0.140	1.83
1	CV (%)	5.18	5.52	6.65

**Recovery:** The recovery of antigen spiked to levels throughout the range of the assay in depleted plasma was evaluated.

Sample	1	2	3	4	
n	4	4	4	4	
Mean (ng/ml)	0.456	0.456 6.58		39.8	
Average %	114	94	107	100	
Recovery					
Pango	113-	92-	98-	97-	
Range	117%	96%	112%	105%	

**Linearity:** To assess the linearity of the assay, samples containing high concentrations of antigen were serially diluted to produce samples with values within the dynamic range of the assay.

Sample	1:2	1:4	1:8	1:16
n	4	4	4	4
Average % of Expected	99	99	106	101
Range	95- 106%	96- 103%	105- 108%	93- 105%

**Sample Values:** Samples were evaluated for the presence of the antigen at varying dilutions.

Sample Type	Dilution	Mean (mg/mL)		
Citrate Plasma	1:1,000,000	10.6		
	1:2,000,000	9.2		

**Specificity:** This assay recognizes total human IgG. Pooled normal plasma from mouse, rat and rabbit was assayed and no significant cross-reactivity was observed. Pooled normal plasma from cyno monkey was assayed and minor cross-reactivity was observed.

# **DISCLAIMER**

This information is believed to be correct but does not claim to be all-inclusive and shall be used only as a guide. The supplier of this kit shall not be held liable for any damage resulting from handling of or contact with the above product.

Example of ELISA Plate Layout
96 Well Plate: 22 Standard wells, 74 Sample wells

_	1	2	3	4	5	6	7	8	9	10	11	12
Α	0	0.2 ng/ml	0.5 ng/ml	1 ng/ml	2 ng/ml	5 ng/ml	10 ng/ml	20 ng/ml	50 ng/ml	100 ng/ml	200 ng/ml	
В	0	0.2 ng/ml	0.5 ng/ml	1 ng/ml	2 ng/ml	5 ng/ml	10 ng/ml	20 ng/ml	50 ng/ml	100 ng/ml	200 ng/ml	
С				<b>0</b> 5								
D			J?									
E												
F												
G												
н												