

# Amplite™ Fluorimetric Goat Anti-Rabbit IgG- HRP Conjugate ELISA Assay Kit \*Red Fluorescence\*

Catalog number: 11541

Unit size: 1000 Tests

Component	Storage	Amount
Component A: Amplite™ Red Peroxidase Substrate	Freeze (<-15 °C), Minimize light exposure	2 vials
Component B: H2O2	Refrigerate (2-8 °C), Minimize light exposure	1 vial (3% stabilized solution, 500 µL)
Component C: Assay Buffer	Freeze (<-15 °C)	1 bottle (100 mL)
Component D: DMSO	Freeze (<-15 °C)	1 vial (1 mL)
Component E: Goat Anti-Rabbit IgG-HRP Conjugate	Freeze (<-15 °C)	1 vial (25 µL)

## OVERVIEW

Horseshoe Peroxidase (HRP) is a small molecule (MW ~40 KD) that is widely used in a variety of biological detections. HRP conjugates are extensively used as secondary detection reagents in ELISAs, immuno-histochemical techniques, Northern, Southern and Western blot analyses. Due to its small size, it rarely causes steric hindrance problem with antibody/antigen complex formation. It is usually conjugated to an antibody in a 4:1 ratio. Additionally, HRP is inexpensive compared to other labeling enzymes. The major disadvantage associated with peroxidase is their low tolerance to many preservatives such as sodium azide that inactivates peroxidase activity even at low concentration. Our Amplite™ Fluorimetric ELISA Assay Kit contains all the essential components including our fluorogenic Amplite™ Red HRP substrate for ELISA detection. The kit provides an optimized assay protocol that is compatible with HTS liquid handling, as little as 10 pg of a Polyclonal antibody in the well of a microplate can be detected. Its signal can be easily read by either fluorescence microplate reader or absorbance microplate reader. It can be used for the assays that detect goat anti-rabbit IgG as the secondary detection agent.

## AT A GLANCE

### Protocol summary

1. Prepare ELISA plate
2. Add Goat Anti-Rabbit IgG-HRP Conjugate working solution (100 µL/well)
3. Incubate at room temperature for 30 minutes
4. Wash the wells (3X with PBS-Tween)
5. Add Peroxidase working solution (100 µL/well)
6. Incubate at room temperature for 30 - 60 minutes
7. Monitor fluorescence intensity at Ex/Em = 540/490 nm (Cutoff = 570 nm)

**Important** Thaw all the kit components at room temperature before starting the experiment.

## KEY PARAMETERS

Instrument: Fluorescence microplate reader  
 Excitation: 540 nm  
 Emission: 590 nm  
 Cutoff: 570 nm  
 Recommended plate: Solid black

Instrument: Absorbance microplate reader  
 Absorbance: 576 ± 5 nm  
 Recommended plate: Clear bottom

## PREPARATION OF STOCK SOLUTIONS

*Unless otherwise noted, all unused stock solutions should be divided into single-use aliquots and stored at -20 °C after preparation. Avoid repeated freeze-thaw cycles.*

### 1. Amplite™ Red Peroxidase Substrate stock solution (200X):

Add 250 µL of DMSO (Component D) into the vial of Amplite™ Red Peroxidase Substrate (Component A) and mix well to make 200X Amplite™ Red Peroxidase Substrate stock solution.

**Note** 50 µL of 200X Amplite™ Red Peroxidase Substrate stock solution is enough for 1 plate. The stock solution should be used promptly. Protect from light.

### 2. H<sub>2</sub>O<sub>2</sub> stock solution (20 mM):

Add 22.7 µL of 3% H<sub>2</sub>O<sub>2</sub> (0.88 M, Component B) into 977µL of Assay Buffer (Component C) to make 20 mM H<sub>2</sub>O<sub>2</sub> stock solution.

**Note** The diluted H<sub>2</sub>O<sub>2</sub> stock solution is not stable. The unused portion should be discarded.

## PREPARATION OF WORKING SOLUTION

Add 50 µL of 200X Amplite™ Red Peroxidase Substrate stock solution and 100 µL of 20 mM H<sub>2</sub>O<sub>2</sub> stock solution into 9.85 mL of Assay Buffer (Component C) and mix well to make Peroxidase working solution. Keep from light.

## PREPARATION OF CELL SAMPLES

For guidelines on cell sample preparation, please visit <https://www.aatbio.com/resources/guides/cell-sample-preparation.html>

## SAMPLE EXPERIMENTAL PROTOCOL

### Prepare ELISA plate:

1. Prepare ELISA microplate (including appropriate controls) by performing all necessary ELISA preparation steps.
2. Add 2 µL of Goat Anti-Rabbit IgG-HRP Conjugate (Component E) into 10 mL of PBS with 1% BSA (PBS-BSA, not included) to make Goat Anti-Rabbit IgG-HRP Conjugate working solution.

**Note** 10 mL of Goat Anti-Rabbit IgG-HRP Conjugate working solution is enough for 1 plate. The concentration of this Goat Anti-Rabbit IgG-HRP Conjugate working solution is recommended as an initial concentration to try. The optimal concentration for each particular application may need to be determined empirically.

3. Wash the ELISA wells three times with PBS containing 0.02% to 0.05% Tween® 20 (PBS-Tween) and drain.

**DISCLAIMER**

AAT Bioquest provides high-quality reagents and materials for research use only. For proper handling of potentially hazardous chemicals, please consult the Safety Data Sheet (SDS) provided for the product. Chemical analysis and/or reverse engineering of any kit or its components is strictly prohibited without written permission from AAT Bioquest. Please call 408-733-1055 or email info@aatbio.com if you have any questions.

4. Add 100  $\mu$ L of Goat Anti-Rabbit IgG-HRP Conjugate working solution into each well.
5. Incubate at room temperature for 30 minutes. Drain off the Goat Anti-Rabbit IgG-HRP Conjugate.
6. Wash the wells three times with PBS-Tween and drain.

**Run Peroxidase assay in ELISA plate:**


1. Add 100  $\mu$ L of Peroxidase working solution into each drained microplate well containing the samples and controls.
2. Incubate the reaction at room temperature for 30 minutes or longer, protected from light.
3. Monitor the fluorescence increase with a fluorescence plate reader at Excitation = 530-570 nm, Emission = 590-600 nm (optimal Ex/Em = 540/590 nm, Cutoff = 570 nm).

**Note** The plate can also be read by an absorbance microplate reader at the wavelength of  $576 \pm 5$  nm. The absorption detection has lower sensitivity compared to fluorescence reading.

**EXAMPLE DATA ANALYSIS AND FIGURES**

The reading (RFU) obtained from the blank standard well is used as a negative control. Subtract this value from the other standards' readings to obtain the baseline corrected values. Then, plot the standards' readings to obtain a standard curve and equation. This equation can be used to calculate Rabbit total IgG samples. We recommend using the Online Four Parameter Logistics Calculator which can be found at:

<https://www.aatbio.com/tools/four-parameter-logistic-4pl-curve-regression-online-calculator>



**Figure 1.** Detection of rabbit total IgG using the Amplite™ Fluorimetric ELISA Kit. Rabbit IgG was diluted into 1  $\mu$ g/mL and made 1 to 3 serial dilutions in 0.2 M sodium bicarbonate buffer at pH 9.4. 100  $\mu$ L/well serial dilutions were coated into a black 96-well plate at 4°C overnight, and blocked with 3% milk in PBS and 0.02% Tween 20 at 4°C overnight. The wells were washed, and assayed using the reagents. 1 to 6000 dilutions of goat anti-rabbit IgG-HRP conjugate were used. The reactions were incubated for 15 to 60 minutes and then measured for fluorescence at Ex/Em = 540/590 nm with Gemini fluorescence microplate reader (Molecular Devices).