

Amplite™ Fluorimetric Beta-Galactosidase Assay Kit *Green Fluorescence*

 Catalog number: 12601
 Unit size: 500 Tests

Component	Storage	Amount
Component A: Fluorescein di-β-D-Galactopyranoside (FDG)	Freeze (< -15 °C), Minimize light exposure	1 vial
Component B: Reaction Buffer	Freeze (< -15 °C), Minimize light exposure	1 bottle (50 mL)
Component C: Stop Buffer	Freeze (< -15 °C), Minimize light exposure	1 vial (25 mL)
Component D: Lysis Buffer	Freeze (< -15 °C), Minimize light exposure	1 vial (25 mL)
Component E: DMSO	Freeze (< -15 °C)	1 vial (500 μL)
Component F: β-Mercaptoethanol	Freeze (< -15 °C), Minimize light exposure	1 vial (500 μL)

OVERVIEW

E. coli beta-galactosidase is a 464 kD tetramer. Each unit of beta-galactosidase consists of five domains, the third of which is the active site. It is an essential enzyme in cells. Deficiencies in this enzyme can result in galactosialidosis or Morquio B syndrome. In E. coli, beta-galactosidase is produced by the activation of LacZ operon. Detection of LacZ expression has become routine to the point of detection of as few as 5 copies of β-galactosidase per cell. This kit uses a fluorogenic galactosidase substrate that can sensitively distinguish LacZ+ vs. LacZ-cells. It can be used either for detecting galactosidase conjugates in ELISA type assay systems or for monitoring LacZ gene expression in cells. The galactosidase-cleaved product has an emission spectra that can be detected with most of fluorescence instruments equipped with a FITC filter set.

AT A GLANCE

Protocol Summary

1. Prepare stable or transient transfected cells with LacZ gene
2. Incubate cells (samples) with test compounds
3. Lyse the cells
4. Transfer the lysate to a microtiter plate
5. Add FDG working solution
6. Incubate at room temperature or 37°C for at least 5 minutes depending on cell type
7. Add stopping solution
8. Monitor fluorescence intensity at Ex/Em = 490/525 nm

Important Thaw all the kit components to room temperature before use.

KEY PARAMETERS

Fluorescence microplate reader

Excitation	490 nm
Emission	525 nm
Cutoff	515 nm
Recommended plate	Solid black

CELL PREPARATION

For guidelines on cell sample preparation, please visit <https://www.aatbio.com/resources/guides/cell-sample-preparation.html>

PREPARATION OF STOCK SOLUTIONS

Unless otherwise noted, all unused stock solutions should be divided into single-use aliquots and stored at -20 °C after preparation. Avoid repeated freeze-thaw cycles.

FDG stock solution (1X)

Add 125 μL of DMSO (Component E) into the vial of FDG (Component A) to make 1X FDG stock solution. Note: 25 μL of FDG is enough for 1 plate. Keep from light.

PREPARATION OF STANDARD SOLUTION

For convenience, use the Serial Dilution Planner:

<https://www.aatbio.com/tools/serial-dilution/12601>

β-Galactosidase standard

Optional (if a standard curve is desired): Prepare a serial dilution of β-galactosidase (E. Coli) standards with 0.3% β- mercaptoethanol assay buffer. Transfer 50 μL aliquot of each point on the standard curve to the control wells of the plate. The highest recommended amount of β-galactosidase is 200 mU/mL (200 - 400 ng). 2X serial dilution of standard curve consisting of 8 points is recommended. Note: Adjust the standard curve to suit the specific experimental conditions, such as cell type, number, transfection efficiency, and size of the culture plates. The dilutions for the standard curve must be prepared freshly each time the assay is performed.

PREPARATION OF WORKING SOLUTION

1. 0.3 % β-mercaptoethanol assay buffer

Add 30 μL of β-mercaptoethanol (Component F) to 10 mL of Reaction Buffer (Component B), and mix well. Note: Additional buffer is needed for preparing enzyme dilution buffer, which is used to generate a standard curve.

2. FDG working solution

Add 25 μL of 1X FDG stock solution into 5 mL of 0.3 % β-mercaptoethanol assay buffer. Note: DO NOT keep FDG solutions at room temperature for an extended period of time as spontaneous hydrolysis will occur.

3. Lysis buffer working solution

Add 5 μL of β-mercaptoethanol (Component F) to 5 mL of Lysis Buffer (Component D) before use. Note: Always add 0.1% β-mercaptoethanol into lysis buffer before lysing the cells

SAMPLE EXPERIMENTAL PROTOCOL

Table 1. Recommended Lysis Buffer working solution volumes for cell culture plates.

Type of culture plates	Lysis Buffer working solutions (μL/well)
96-well plate	50
24-well plate	250
12-well plate	500
6-well plate	1000
60 mm plate	2000
100 mm plate	4000

Prepare cell extracts from mammalian cells

1. Treat cells containing LacZ gene with test compounds for a desired period of time.
2. Wash the cells twice with 1X PBS. Do not dislodge the cells.

- Lyse cells accordingly with Lysis Buffer working solution. For adherent cells: Add Lysis Buffer working solution to the culture plates. See table 1 for recommended volumes. For non-adherent cells: Pellet the cells into centrifuge tube, and add 50 - 2000 μL (depending on the size of the cell pellet) of Lysis Buffer working solution to the tube.
- Incubate cells from previous step at room temperature for 10 - 15 minutes, and gently swirl the plates or tubes several times to ensure complete lysis.
- Proceed to the FDG assay or freeze the sample at $-80\text{ }^{\circ}\text{C}$ until use. *Note:* A good lysis can also be obtained by a quick freeze-and-thaw cycle (freeze 1 - 2 hours at $-20\text{ }^{\circ}\text{C}$ to $-80\text{ }^{\circ}\text{C}$ and thaw at room temperature). Alternatively, centrifuge the cell lysis for 2 - 3 minutes to pellet the insoluble material, and then assay the supernatant.

Run β -galactosidase assay

- Thaw the tube or plate of lysed cells at room temperature if needed. Perform the assay directly on the 96-well plate if the cells were seeded in a 96-well plate.
- Add 50 μL of cell extracts into each well of the 96-well plate. Save some control wells for the standard curve (50 μL /well) if a standard curve is desired. *Note:* If necessary, dilute the lysate in Lysis Buffer working solution when transfection efficiency is very high or reduce the volume of lysis buffer when transfection efficiency is low. If the transfection is performed in a 96-well plate, or a stable cell line was seeded into a 96-well plate, perform the assay directly on the plate. For endogenous β -galactosidase activity control, add 50 μL of cell lysate from non-transfected cells. For blank control, add 50 μL of 1X lysis buffer.
- Add 50 μL of FDG working solution to each well. Incubate the plate at room temperature or $37\text{ }^{\circ}\text{C}$ for approximately 5 min to 4 hr depending on the cell type.
- Add 50 μL of Stop Buffer (Component C) to each well. The stop buffer causes an increase in the fluorescence intensity of the product, in addition to terminate the reaction.
- Measure the fluorescence intensity of the solution in each well with a fluorescence microplate reader at $\text{Ex/Em} = 490/525\text{ nm}$.
- Quantify β -galactosidase expression based on a linear standard curve.

EXAMPLE DATA ANALYSIS AND FIGURES




Figure 1. β -galactosidase dose response was measured with Amplitude™ Fluorimetric beta-Galactosidase Assay Kit in a Costar solid black 96-well plate using Gemini fluorescence microplate reader (Molecular Devices).

DISCLAIMER

AAT Bioquest provides high-quality reagents and materials for research use only. For proper handling of potentially hazardous chemicals, please consult the Safety Data Sheet (SDS) provided for the product. Chemical analysis and/or reverse engineering of any kit or its components is strictly prohibited without written permission from AAT Bioquest. Please call 408-733-1055 or email info@aatbio.com if you have any questions.