

Amplite™ Fluorimetric cADP-Ribose Assay Kit

Catalog number: 20305
Unit size: 100 tests

| Component | Storage | Amount |
|-------------------------------------|---|-------------------|
| Component A: Quest Fluor™ NAD Probe | Freeze (<-15 °C), Minimize light exposure | 1 bottle (5 mL) |
| Component B: ADPRC Enzyme mix | Freeze (<-15 °C), Minimize light exposure | 1 vial |
| Component C: Assay Solution I | Freeze (<-15 °C) | 1 bottle (12 mL) |
| Component D: Assay Solution II | Freeze (<-15 °C) | 1 bottle (5 mL) |
| Component E: Enhancer Solution | Freeze (<-15 °C), Minimize light exposure | 1 bottle (3.5 mL) |
| Component F: cADPR Standard | Freeze (<-15 °C), Minimize light exposure | 1 vial |

OVERVIEW

cADP-ribose (cADPR) is a novel Ca²⁺ messenger derived from NAD⁺. ADP-riboxyl cyclase (ADPRC) catalyzes the synthesis of cADPR from NAD⁺, but the reaction can be reversed in the presence of high concentration of nicotinamide, producing NAD⁺ from cADPR stoichiometrically. The resultant NAD⁺ can be detected using our newly developed NAD sensor Quest Fluor™ NAD reagent. This makes monitoring cADPR in tissues and cell cultures possible in the low nM range. The NAD⁺ detection using Quest Fluor™ NAD reagent is specific to NAD⁺ and has no reaction to NADH. The fluorescent signal can be readily detected. This assay can be performed in a convenient 96-well or 384-well microtiter plate.

AT A GLANCE

Protocol summary

1. Prepare and add cADPR standards and/or test samples (50 µL)
2. Prepare and add ADPRC working solution (50 µL)
3. Incubate at room temperature for 1 hour
4. Add 40 µL Quest Fluor™ NAD Probe
5. Add 40 µL Assay Solution
6. Incubate at room temperature for 20 minutes
7. Add 30 µL Enhancer Solution
8. Incubate at room temperature for 20 minutes
9. Monitor fluorescence intensity at Ex/Em = 420/480 nm

Important Thaw all the kit components at room temperature before starting the experiment.

KEY PARAMETERS

| | |
|--------------------|--------------------------------|
| Instrument: | Fluorescence microplate reader |
| Excitation: | 420 nm |
| Emission: | 480 nm |
| Cutoff: | 435 nm |
| Recommended plate: | Solid black |

PREPARATION OF STOCK SOLUTIONS

Unless otherwise noted, all unused stock solutions should be divided into single-use aliquots and stored at -20 °C after preparation. Avoid repeated freeze-thaw cycles.

cADPR standard stock solution (5 mM):

Add 10 µL of ddH₂O into the vial of cADPR standard (Component F) and mix them well.

Note The unused cADPR stock solution should be stored at -20°C in single use aliquots.

PREPARATION OF STANDARD SOLUTION

cADPR standard

For convenience, use the Serial Dilution Planner:

<https://www.aatbio.com/tools/serial-dilution/20305>

Note Prepare cADPR serial dilutions in Assay Solution I (Component C).

PREPARATION OF WORKING SOLUTION

ADPRC working solution:

Add 50 µL ddH₂O into the vial of ADPRC Enzyme Mix (Component B) and mix well. Transfer whole content into 5 mL of Assay Solution I (Component C) and mix them well.

Note ADPRC working solution is not stable, use it promptly.

SAMPLE EXPERIMENTAL PROTOCOL

Table 1. Layout of cADPR standards and test samples in a solid black 96-well microplate. ST = cADPR standard (ST1-ST7); BL = blank control; TS = test sample.

| | | | |
|-----|-----|-----|-----|
| BL | BL | TS | TS |
| ST1 | ST1 | ... | ... |
| ST2 | ST2 | ... | ... |
| ST3 | ST3 | | |
| ST4 | ST4 | | |
| ST5 | ST5 | | |
| ST6 | ST6 | | |
| ST7 | ST7 | | |

Table 2. Reagent composition for each well.

| Well | Volume | Reagent |
|---------|--------|--------------------------------|
| ST1-ST7 | 50 µL | Serial Dilution |
| BL | 50 µL | Assay Solution I (Component C) |
| TS | 50 µL | Test Sample |

NAD generation assay

1. Add 50 µL of cADPR standard, blank control, and test samples to solid black 96-well microplate (As shown in Table 1 and Table 2).
2. Add 50 µL/well of ADPRC working solution into each well of cADPR standard, blank control and test samples. *Note:* For a 384-well plate, add 12.5 µL of sample and 12.5 µL of ADPRC Reaction Mix Solution into each well.

3. Incubate the reaction at room temperature for 60 minutes, protected from light.

NAD detection assay

1. Add 40 μL Quest Fluor™ NAD Probe (Component A) into each well of cADPR standard, blank control, and test samples (total of 140 $\mu\text{L}/\text{well}$), mix well.
2. Add 40 μL Assay Solution II (Component D) into each well (total of 180 $\mu\text{L}/\text{well}$), mix well.

Note For a 384-well plate, add 10 μL of Quest Fluor™ NAD Probe and 10 μL Assay Solution II into each well.

3. Incubate the reaction at room temperature for 10 - 20 minutes, protected from light.
4. Add 30 μL Enhancer Solution (Component E) to each well to make the total NAD assay volume of 210 $\mu\text{L}/\text{well}$, and incubate at room temperature for 10-20 minutes, protected from light.

Note For a 384-well plate, add 7.5 μL Enhancer Solution.

5. Monitor the fluorescence increase with a fluorescence plate reader at 420/480 nm.

EXAMPLE DATA ANALYSIS AND FIGURES

The reading (RFU) obtained from the blank standard well is used as a negative control. Subtract this value from the other standards' readings to obtain the baseline corrected values. Then, plot the standards' readings to obtain a standard curve and equation. This equation can be used to calculate cADPR samples. We recommend using the Online Linear Regression Calculator which can be found at:

<https://www.aatbio.com/tools/linear-logarithmic-semi-log-regression-online-calculator>




Figure 1. The concentration of cADPR was detected using Amplite™ Fluorimetric cADP-Ribose Assay Kit. Different concentrations of cADPR were incubated with ADPRC reaction mix for 20 min at RT before NAD detection reagent was added. The lowest detected concentration of cADPR is 100 nM.

DISCLAIMER

AAT Bioquest provides high-quality reagents and materials for research use only. For proper handling of potentially hazardous chemicals, please consult the Safety Data Sheet (SDS) provided for the product. Chemical analysis and/or reverse engineering of any kit or its components is strictly prohibited without written permission from AAT Bioquest. Please call 408-733-1055 or email info@aatbio.com if you have any questions.