


Modification of Peptides, Oligonucleotides and Other Small Biomolecules with Tide Fluor™ (TF) Dye Maleimides

Introduction

The sulfonated cyanine dyes have improved fluorescence quantum yield than those of non-sulfonate cyanines. They are commercially available as Cy3®, Cy5®, Cy5.5® and Cy7®. However, these commercial protein-labeling CyDyes and Alexa Fluor® dyes are cost-prohibitive for labeling peptides and oligonucleotides that require much larger amount of labeling dyes. The improved and affordable Tide Fluor™ dyes are optimized as building blocks for developing FRET oligonucleotides and peptides. Our Tide Fluor™ dyes have stronger fluorescence and higher photostability than the classic fluorophores (such as fluoresceins, rhodamines and cyanines). They are the best affordable fluorescent dyes for labeling peptides and oligonucleotides without compromised performance.

- ➔ Optimized to label peptides and nucleotides with performance comparable to Alexa Fluor® Dyes
- ➔ Optimized to pair with Tide Quencher™ dark acceptors to maximize FRET efficiency
- ➔ pH-insensitive and environment-insensitive fluorescence for developing robust assays
- ➔ Higher photostability to improve the quality of fluorescence imaging
- ➔ A variety of reactive forms available for conjugations

Labeling Mechanism



Because free thiol (SH) groups, also called mercapto groups, are not present as abundantly as amino groups in biological samples, thiol-reactive reagents often provide a means of selectively modifying a biomolecule at a defined site. Thiol-reactive dyes are often used to prepare fluorescent peptides and oligonucleotides for probing biological structures, functions and interactions. There are many types of thiol-reactive dyes reported in the literature, including iodoacetamides, disulfides, maleimides, vinyl sulfones and various electron-deficient aryl halides and sulfonates. Maleimide is by far the most popular thiol-reactive moiety. Maleimides readily react with thiol moieties of biopolymers to form thioether conjugates even under neutral conditions. The thioether bond formed is quite stable. Maleimides require conjugation conditions less stringent than those of iodoacetamides, and do not react with histidine and methionine under physiological conditions. For example, most conjugations can be done at room temperature at neutral pH.

Storage and Handling

The dye labeled oligos should be stored at $\leq -15^{\circ}\text{C}$, and kept from light. For longer storage, dye labeled oligo could be divided as single-used aliquot and stored at $\leq -15^{\circ}\text{C}$. PROTECT FROM LIGHT.

Sample Protocols (FOR REFERENCE PURPOSE ONLY)

Always wear the protection apparatus (such as gloves). Warm all the components before opening them, and immediately prepare the required solutions before starting the conjugation. Avoid repeated freezing and thawing if possible. Any solutions containing the dye should be kept from light.

Label Thiol-Modified Oligonucleotides with Tide Fluor™ Dyes

The following protocol has been optimized for labeling 200 µg (~6 $A_{260\text{ nm}}$ units) of a proprietary oligonucleotide. You need modify the protocol to get the best results for your particular application by multiple experimentations. **YOUR THIOL-MODIFIED OLIGO MUST BE TREATED TO REMOVE SMALL THIOL COMPOUNDS (SUCH AS DTT) THAT RAPIDLY REACTS AND CONSUMES DYE MALEIMIDES.**

1. Prepare Oligo Solution (Solution A)

Dissolve your thiol-modified oligo (~200 µg) in a PBS buffer (100 µL, pH 7.2).

Note 1: The oligonucleotide must be synthesized with a thiol group. See Appendix 1 for the purification of thiol-modified oligos.

Note 2: Avoid buffers that contain thiol compounds, such as DTT, as these compete for conjugation with the thiol-reactive compound.

2. Prepare Dye Solution (Solution B)

Dissolve 1 mg dye maleimide in 100 µL DMSO (>10 mg/mL if possible) by pipetting up and down. Centrifuge the solution stock on the sides of the vial to the vial bottom.

Note: prepare the DMSO dye solution before starting the conjugation. Extended storage of the dye solution may reduce the dye activity. Any solutions containing the dye should be kept from light.

3. Run Conjugation Reaction

3.1 To the dye solution (B, 20-50 µL) add the oligo solution (A, 100 µL) with stirring or shaking (keeping the reaction mixture from light).

3.2 Rotate or shake the reaction mixture for 4-6 hours at room temperature on a rotator or shaker.

Note: Gently vortex tap the vial every 10 minutes for the first hour to ensure that the reaction solution remains well mixed. Do not mix violently, as material may be left on the sides of the vial. After six hours, 50–90% of the thiol-modified oligonucleotide molecules should be labeled. The reaction might be incubated overnight if it is more convenient. However, overnight incubation will not result in a greater labeling efficiency in most cases.

4. Purify Dye-Oligo Conjugate

4.1 Preliminary purification by ethanol precipitation of labeled oligonucleotide

a. Add 20 µL (one-tenth reaction solution volume in general) of 3 M NaCl and 300 µL cold absolute ethanol (two and half reaction solution volume volumes in general) to the reaction vial.

b. Mix the solution well and place it at –20°C for 30 minutes.

c. Centrifuge the solution in a microcentrifuge at 10,000 to 15,000 × g for 30 minutes.

Note: Loss of sample may occur if the centrifugation is not long enough.

d. Carefully remove the supernatant, rinse the pellet 1-3 times with cold 70% ethanol and dry briefly.

Note: Some unreacted labeling reagent may have precipitated over the course of the reaction or may be stuck on the walls of the reaction vial. This material should be completely redissolved by extensive vortex mixing before centrifugation. Redissolving the labeling reagent ensures that the precipitated oligonucleotide will be minimally contaminated with unreacted label.

4.2 Final Purification by HPLC or by gel electrophoresis

See Appendix I

Label Peptides with Tide Fluor™ Dyes

The following protocol has been optimized for labeling 10 mg of a proprietary peptide (MW ~ 2000) that contains only a single free thiol group. **YOU NEED MODIFY THE PROTOCOL TO ARCHIE THE BEST RESULTS FOR YOUR PARTICULAR APPLICATION BY MULTIPLE EXPERIMENTATIONS.**

1. **Prepare Peptide Solution (Solution A)**

Dissolve the peptide (to be labeled) in DMF, DMSO or water (>10 mg/mL if possible).

Note: For effective labeling, higher peptide concentration is preferable.

2. **Prepare Dye Solution (Solution B)**

Dissolve the dye maleimide in DMF or DMSO (>5 mg/mL if possible).

Note: For effective labeling, higher dye concentration is preferable.

3. **Mix the peptide and dye solutions (from Steps 1 and 2) in roughly equal molar ratio.**

Note: Either peptide or dye maleimide can be used in excess depending on the cost of peptide and dye. Good separation of free dye, free peptide and dye-peptide conjugate is another important factor that need be considered as well to decide that either peptide or dye need be used in slight excess.

4. **Stir the reaction mixture at room temperature for 12-24 hours.**

Note: This reaction is quite efficient. It is generally completed within a few hours. If the reaction does not generate a considerable amount of desired conjugate within 2 hours, you need check the reaction pH, a very critical factor. The optimal pH is 4-6. Too low pH significantly decreases the reaction rate while in high pH solutions peptides tend to be oxidized to give the disulfides. Peptides often contain TFA, therefore it is important to check the reaction mixture pH and adjust it to pH 4-6 using NaHCO₃ if necessary.

5. **Purify Dye-Peptide Conjugate**

The reaction solution was concentrated and purified on a C18 column to afford the desired conjugate. The fractions were analyzed by HPLC, and the fractions of >97% purity were pooled and lyophilized.

Note 1: HPLC Purification Conditions: TEAB buffer (triethyl ammonium bicarbonate, 0.25 mmol, pH=7.0-8.0) was used as buffer A and acetonitrile as buffer B. The HPLC was run from 0% B to 30% B in 60 min (flow rate: 100 mL/min).

Note 2: Avoid strong light during the operation.

References

1. Hermanson GT (1996). *Biocojugate Techniques*, Academic Press, New York.
2. Sambrook J., Fritsch E.F. and Maniatis, T. (1989). *Molecular Cloning: A Laboratory Manual, Second Edition*, Cold Spring Harbor Laboratory.
3. Brinkley M (1992). A brief survey of methods for preparing protein conjugates with dyes, haptens, and cross-linking reagents. *Bioconjug Chem* **3**, 2-13.

Appendix I. Chemical Properties of Tide Fluor™ Fluorescent Labeling Dyes

Tide Fluor™ dyes have improved labeling performance than the classic fluorescent labeling dyes such as FITC, TRITC, Texas Red®, Cy3, Cy5 and Cy7. They are the best affordable fluorescent dyes (alternative to Alexa Fluor® dyes) for labeling oligos and peptides without compromised performance.

Labeling Dye	Cat#	Product Description	Reactivity	Adduct MW Calculation*
TF1	2236	Tide Fluor™ 1 azide [TF1 azide]	Azide	+ 301
	2237	Tide Fluor™ 1 alkyne [TF1 alkyne]	Alkyne	+ 270
	2238	Tide Fluor™ 1 acid [TF1 acid]	NH ₂ and OH	+ 215
	2239	Tide Fluor™ 1 amine [TF1 amine]	CO ₂ H	+ 257
	2242	Tide Fluor™ 1 maleimide [TF1 maleimide]	SH	+ 355
	2244	Tide Fluor™ 1 succinimidyl ester [TF1 SE]	Aliphatic amine	+ 215
TF2	2245	Tide Fluor™ 2 acid [TF2 acid]	NH ₂ and OH	+ 469
	2246	Tide Fluor™ 2 amine [TF2 amine]	CO ₂ H	+ 398
	2247	Tide Fluor™ 2 maleimide [TF2 maleimide]	SH	+ 680
	2248	Tide Fluor™ 2 succinimidyl ester [TF2 SE]	Aliphatic amine	+ 469
	2252	Tide Fluor™ 2 azide [TF2 azide]	Azide	+ 555
	2253	Tide Fluor™ 2 alkyne [TF2 alkyne]	Alkyne	+ 524
TF2WS	2348	Tide Fluor™ 2WS acid [TF2WS acid]	NH ₂ and OH	+ 1095
	2349	Tide Fluor™ 2WS succinimidyl ester [TF2WS SE]	Aliphatic amine	+ 1095
TF3	2254	Tide Fluor™ 3 azide [TF3 azide]	Azide	+ 526
	2255	Tide Fluor™ 3 alkyne [TF3 alkyne]	Alkyne	+ 495
	2268	Tide Fluor™ 3 acid [TF3 acid]	NH ₂ and OH	+ 440
	2269	Tide Fluor™ 3 amine [TF3 amine]	CO ₂ H	+ 496
	2270	Tide Fluor™ 3 maleimide [TF3 maleimide]	SH	+ 580
	2271	Tide Fluor™ 3 succinimidyl ester [TF3 SE]	Aliphatic amine	+ 440
TF3WS	2344	Tide Fluor™ 3WS maleimide [TF3WS maleimide]	SH	+ 846
	2345	Tide Fluor™ 3WS acid [TF3WS acid]	NH ₂ and OH	+ 706
	2346	Tide Fluor™ 3WS succinimidyl ester [TF3WS SE]	Aliphatic amine	+ 706
TF4	2285	Tide Fluor™ 4 acid [TF4 acid]	NH ₂ and OH	+ 544
	2286	Tide Fluor™ 4 amine [TF4 amine]	CO ₂ H	+ 586
	2287	Tide Fluor™ 4 maleimide [TF4 maleimide]	SH	+ 755
	2289	Tide Fluor™ 4 succinimidyl ester [TF4 SE]	Aliphatic amine	+ 544
	2300	Tide Fluor™ 4 azide [TF4 azide]	Azide	+ 630
	2301	Tide Fluor™ 4 alkyne [TF4 alkyne]	Alkyne	+ 599
TF5WS	2275	Tide Fluor™ 5WS azide [TF5WS azide]	Azide	+ 1078
	2276	Tide Fluor™ 5WS alkyne [TF5WS alkyne]	Alkyne	+ 787
	2278	Tide Fluor™ 5WS acid [TF5WS acid]	NH ₂ and OH	+ 732
	2279	Tide Fluor™ 5WS amine [TF5WS amine]	CO ₂ H	+ 774
	2280	Tide Fluor™ 5WS maleimide [TF5WS maleimide]	SH	+ 873
	2281	Tide Fluor™ 5WS succinimidyl ester [TF5WS SE]	Aliphatic amine	+ 732

[Continued on next page]

TF6WS	2291	Tide Fluor™ 6WS acid [TF6WS acid]	NH ₂ and OH	+ 899
	2292	Tide Fluor™ 6WS amine [TF6WS amine]	CO ₂ H	+ 941
	2293	Tide Fluor™ 6WS maleimide [TF6WS maleimide]	SH	+ 1039
	2294	Tide Fluor™ 6WS succinimidyl ester [TF6WS SE]	Aliphatic amine	+ 899
	2302	Tide Fluor™ 6WS azide [TF6WS azide]	Azide	+ 1079
	2303	Tide Fluor™ 6WS alkyne [TF6WS alkyne]	Alkyne	+ 1048
TF7WS	2304	Tide Fluor™ 7WS azide [TF7WS azide]	Azide	+ 845
	2305	Tide Fluor™ 7WS alkyne [TF7WS alkyne]	Alkyne	+ 813
	2330	Tide Fluor™ 7WS acid [TF7WS acid]	NH ₂ and OH	+ 758
	2331	Tide Fluor™ 7WS amine [TF7WS amine]	CO ₂ H	+ 801
	2332	Tide Fluor™ 7WS maleimide [TF7WS maleimide]	SH	+ 899
	2333	Tide Fluor™ 7WS succinimidyl ester [TF7WS SE]	Aliphatic amine	+ 758
TF8WS	2306	Tide Fluor™ 8WS azide [TF8WS azide]	Azide	+ 1011
	2307	Tide Fluor™ 8WS alkyne [TF8WS alkyne]	Alkyne	+ 980
	2335	Tide Fluor™ 8WS acid [TF8WS acid]	NH ₂ and OH	+ 925
	2336	Tide Fluor™ 8WS amine [TF8WS amine]	CO ₂ H	+ 967
	2337	Tide Fluor™ 8WS maleimide [TF8WS maleimide]	SH	+ 1065
	2338	Tide Fluor™ 8WS succinimidyl ester [TF8WS SE]	Aliphatic amine	+ 925

* The molecular weight of the desired conjugate = the molecular weight of the free unlabeled molecule + the value listed in the table.

Appendix II. Spectral Properties of Tide Fluor™ Fluorescent Labeling Dyes

Labeling Dye	Extinction Coefficient¹ (cm⁻¹M⁻¹)	Abs (nm)	Em (nm)	FQY²	CF at 260 nm³	CF at 280 nm⁴
TF1	20,000	345	442	0.95	0.246	0.187
TF2	75,000	500	527	0.90	0.288	0.201
TF2WS	75,000	502	525	0.90	0.211	0.091
TF3	85,000	555	584	0.85	0.331	0.201
TF3WS	150,000	555	565	0.10 ⁵	0.079	0.079
TF4	90,000	590	618	0.91	0.489	0.436
TF5WS	250,000	649	664	0.25	0.023	0.027
TF6WS	220,000	676	695	0.18	0.111	0.009
TF7WS	275,000	749	775	0.12	0.009	0.049
TF8WS	250,000	775	807	0.08	0.103	0.109

Notes: 1. Extinction Coefficient at their maximum absorption wavelength; 2. FQY = fluorescence quantum yield in aqueous buffer (pH 7.2); 3. CF at 260 nm is the correction factor used for eliminating the dye contribution to the absorbance at 260 nm (for oligo and nucleic acid labeling); 3. CF at 280 nm is the correction factor used for eliminating the dye contribution to the absorbance at 280 nm (for peptide and protein labeling); 5. Fluorescence intensity is significantly increased upon coupled to proteins or long peptides.

Appendix III. FRET Selection Guide of Tide Fluor™ Dyes

Tide Fluor™ Donor	Ex(nm)	Em (nm)	Features and Benefits	Ordering Information
Tide Fluor™ 1 (TF1)	345 nm	442 nm	<i>Alternative to EDANS</i> <ul style="list-style-type: none"> • Much stronger absorption • Much stronger fluorescence • Less environment-sensitive 	#2236 & 2237 (TF1 Click chemistry) #2238 (TF1 acid) #2239 (TF1 amine) #2242 (TF1 maleimide, SH-reactive) #2244 (TF1 SE, NH ₂ -reactive)
Tide Fluor™ 2 (TF2) Tide Fluor™ 2WS (TF2WS)	500 nm 502 nm	527 nm 525 nm	<i>Alternative to FAM, FITC and Alexa Fluor® 488</i> pH-insensitive fluorescence Photostable	#2245 (TF2 acid) & 2348 (TF2WS acid) #2246 (TF2 amine) #2247 (TF2 maleimide, SH-reactive) #2248 (TF2, SE) & #2249 (TF2WS SE) #2252 & 2253 (Click chemistry)
Tide Fluor™ 3 (TF3) Tide Fluor™ 3WS (TF3WS)	555 nm 555 nm	584 nm 565 nm	<i>Alternative to Cy3® and Alexa Fluor® 555</i> Strong fluorescence Photostable	#2254 & 2255 (TF3 Click chemistry) #2268 (TF3 acid) & 2345 (TF3WS acid) #2269 (TF3 amine) #2270 (TF3 maleimide, SH-reactive) #2271 (TF3 SE) & #2346 (TF3WS SE)
Tide Fluor™ 4 (TF4)	590 nm	618 nm	<i>Alternative to ROX, Texas Red® and Alexa Fluor® 594</i> Strong fluorescence Photostable	#2285 (TF4 acid) #2286 (TF4 amine) #2287 (TF4 maleimide, SH-reactive) #2289 (TF4 SE, NH ₂ -reactive) #2300 & 2301 (TF4 Click chemistry)
Tide Fluor™ 5WS (TF5WS)	649 nm	664 nm	<i>Alternative to Cy5® and Alexa Fluor® 647</i> Strong fluorescence Photostable	#2275 & 2276 (TF5WS Click chemistry) #2278 (TF5WS acid) #2279 (TF5WS amine) #2280 (TF5WS maleimide, SH-reactive) #2281 (TF5WS SE, NH ₂ -reactive)
Tide Fluor™ 6WS (TF6WS)	676 nm	695 nm	<i>Alternative to Cy5.5®, IRDye® 700 and Alexa Fluor® 680</i> Strong fluorescence Photostable	#2291 (TF6WS acid) #2292 (TF6WS amine) #2293 (TF6WS maleimide, SH-reactive) #2294 (TF6WS SE, NH ₂ -reactive) #2302 & 2303 (TF6WS Click chemistry)
Tide Fluor™ 7WS (TF7WS)	749 nm	775 nm	<i>Alternative to Cy7® and Alexa Fluor® 750</i> Strong fluorescence Photostable	#2304 & 2305 (TF7WS Click chemistry) #2330 (TF7WS acid) #2331 (TF7WS amine) #2332 (TF7WS maleimide, SH-reactive) #2333 (TF7WS SE, NH ₂ -reactive)
Tide Fluor™ 8WS (TF8WS)	775 nm	807 nm	<i>Alternative to IRDye® 800</i> Stronger fluorescence Higher Photostability	#2306 & 2307 (TF4 Click chemistry) #2335 (TF8WS acid) #2336 (TF8WS amine) #2337 (TF8WS maleimide, SH-reactive) #2338 (TF8WS SE, NH ₂ -reactive)

*Texas Red® and Alexa Fluor® are the trademarks of Molecular Probes. CyDye, Cy3®, Cy5®, Cy5.5® and Cy7® are the trademarks of GE Health Care. IRDye® 700 and IRDye® 800 are the trademarks of Li-COR. Tide Fluor™ is the trademark of AAT Bioquest.

Appendix IV. HPLC Purification of Dye Oligonucleotide Conjugates

Ethanol Precipitation

Some commercial oligonucleotides often contain some interfering compounds, especially amines, such as triethylamine or Tris, and ammonium salts, we strongly recommend you to extract and precipitate the commercial oligo samples prior to initiating your labeling reaction. On the other hand, the labeling mixture contains labeled oligonucleotide, unlabeled oligonucleotide, hydrolyzed dye acid and unincorporated dye SE. The impurities of hydrolyzed dye acid and unincorporated dye SE resulted from the labeling reaction can be effectively removed by ethanol precipitation. The following protocol was optimized for the further purification of 0.1–1 mg commercial oligonucleotide sample that was purified by HPLC (3–30 A260 units).

- 1) Dissolve your target oligonucleotide in 100 μL of deionized water and extract three times with an equal volume of chloroform.
- 2) Precipitate the oligonucleotide by adding one-tenth volume (10 μL) of 3 M NaCl and two and a half volumes (250 μL) of cold absolute ethanol. Mix well and place at -20°C for 30 minutes.
- 3) Centrifuge the solution in a microcentrifuge at 10,000 to 15,000 g for 30 minutes.
- 4) Carefully remove the supernatant, rinse the pellet 1-3 times with cold 70% ethanol, and dry under a vacuum.
- 5) Dissolve the dry pellet in deionized water to achieve a final concentration of $>50 \mu\text{g}/\mu\text{L}$. This amine-modified oligonucleotide stock solution may be immediately used or stored frozen at $\leq -15^{\circ}\text{C}$.

Purification by HPLC

Labeled oligonucleotides can be purified by reverse-phase HPLC using a standard analytical C8 or C18 column using an analytical or semi-preparative HPLC instrument. The following protocol was optimized for the further purification of 0.1–1 mg labeled oligonucleotide (3–30 A260 units).

- 1) Dissolve the pellet from the ethanol precipitation in 0.1 M triethylammonium acetate (TEAA).
- 2) Load the dissolved pellet onto the column in 0.1 M TEAA and run a linear 5–95% acetonitrile gradient over 30 minutes.
Note 1: There will be peaks that correspond to the unlabeled oligonucleotide, the labeled oligonucleotide, and the free dye. The actual order and number of these peaks depends on the length of the oligonucleotide and the purity of the sample.
Note 2: To determine the identity of the peaks, monitor the absorbance at both 260 nm and at the absorbance maxima (λ_{max}) for the dye. For instruments with only one detector, two small samples should be run, each monitored at a different wavelength. Unlabeled oligonucleotide will show an absorbance at 260 nm only. Both the free dye and the labeled oligonucleotide will have absorbance at both 260 nm (A260 for oligo) and at the absorbance maximum of the dye (Amax for dye); The dye-labeled oligonucleotide will have a higher A260:Amax ratio than the dye or hydrolyzed dye.

Purification by Gel Electrophoresis

- 1) Pour a 0.5 mm-thick polyacrylamide slab gel.
Note: For oligonucleotides less than 25 bases in length, use 19% acrylamide, for oligonucleotides 25–40 bases, 15% acrylamide, and for oligonucleotides 40–100 bases, 12% acrylamide.
- 2) Resuspend the pellet from ethanol precipitation in 200 μL of 50% formamide, and incubate at 55°C for 5 minutes to disrupt any secondary structure.
- 3) Load the warmed oligonucleotide onto the gel and load an adjacent well with 50% formamide plus 0.05% bromophenol blue. The bromophenol blue will migrate at approximately the same rate as the oligonucleotide.
Note: You may need to use several wells.
- 4) Run the gel until the bromophenol blue indicator dye is two-thirds of the way down the gel.
- 5) Remove the gel from the glass plates and place on Saran Wrap.
- 6) Lay the gel on a fluorescent TLC plate.
- 7) Locate the labeled and unlabeled oligonucleotides by illumination with a handheld UV source.