

Labeling Alkyne-Modified Biomolecules with Fluorescent Dye Azides

Labeling Oligonucleotides with Dye Azides

1. Prepare the following stock solutions:
 - 200 mM THPTA [tris(3-hydroxypropyltriazolylmethyl)amine] in water
 - 100 mM CuSO₄ in water
 - Alkyne-modified oligo in water (as concentrated as possible, e.g., >10 mg/mL)
 - 100 mM sodium ascorbate in water
 - 10 mM dye azide in DMSO or water (see our website for recommended solvent)
2. Mix and vortex well CuSO₄ with THPTA in a 1:2 ratio for several minutes before the reaction. This working solution is stable for several weeks when frozen.
3. To the alkyne-modified oligo solution, add an excess of dye azide (2-5 equivalents by molar ratio).
4. Add 5 equivalents of THPTA/CuSO₄ working solution (from Step 1)
5. Add 10-30 equivalents of sodium ascorbate.
6. Stir, vortex or shake the reaction mixture at room temperature for 30-60 minutes.
7. Ethanol-precipitate or purify the oligo by your desired method (e.g., HPLC).

Labeling Peptides with Dye Azides

1. Prepare the following stock solutions:
 - 200 mM THPTA ligand in water
 - 100 mM CuSO₄ in water
 - Alkyne-modified peptide in water or DMF (depending on your peptide solubility, >10 mg/mL if possible)
 - 100 mM sodium ascorbate in water
 - 10 mM dye azide in DMSO or water (see our website for recommended solvent)
2. Incubate CuSO₄ with THPTA ligand in a 1:2 ratio several minutes before the reaction. This solution is stable for several weeks when frozen.
3. To the alkyne-modified peptide solution, add an excess of dye azide (5-10 equivalents by molar ratio).
4. Add 5-10 equivalents of THPTA/CuSO₄.
5. Add 10-20 equivalents of sodium ascorbate.
6. Stir, vortex or shake the reaction mixture at room temperature for 30-60 minutes.
7. Purify your desired peptide by HPLC.

Labeling Small Organic Alkyne Molecules with Dye Azides

1. Prepare the following stock solutions:
 - 200 mM THPTA ligand in water
 - 100 mM CuSO₄ in water
 - Alkyne compound in water or DMF (depending on your compound solubility, >10 mg/mL if possible,)
 - 100 mM sodium ascorbate in water
 - 10 mM dye azide in DMSO or water (see our website for recommended solvent).
2. Incubate CuSO₄ with THPTA ligand in a 1:2 ratio several minutes before the reaction. This solution is stable for several weeks when frozen.
3. To the alkyne solution, add an excess of dye azide (5-10 equivalents by molar ratio).
4. Add 25 equivalents of THPTA/CuSO₄.
5. Add 50 equivalents of sodium ascorbate.
6. Stir the reaction mixture at room temperature for 30-60 minutes.
7. Purify your desired molecule by chromatography or other methods.

Labeling Biopolymers with Dye Azides

1. Prepare the following stock solutions:
 - 200 mM THPTA ligand in water
 - 100 mM CuSO₄ in water
 - Alkyne-modified biopolymer in water (as concentrated as possible, e.g., >5 mg/mL)
 - 100 mM sodium ascorbate in water
 - 10 mM dye azide in DMSO or water (see our website for recommended solvent).
2. Incubate CuSO₄ with THPTA ligand in a 1:2 ratio several minutes before the reaction. This solution is stable for several weeks when frozen.
3. To the alkyne-modified biopolymer solution, add an excess of dye azide (Loading ratio: 5-20 dye azide/alkyne).
4. Add 5 molar equivalents (referenced to dye azide) of THPTA/CuSO₄.
5. Add 10 equivalents of sodium ascorbate (referenced to dye azide).
6. Stir, vortex or shake the reaction mixture at room temperature for 30-60 minutes.
7. Purify your desired molecule by gel filtration or dialysis.

Labeling Cells, Cell Lysates or Biological Samples with Dye Azides or Dye Alkynes

1. Prepare the following click solutions:
 - 100 mM THPTA ligand in aqueous buffer or water
 - 20 mM CuSO₄ in water
 - 300 mM sodium ascorbate in water
 - 2.5 mM alkyne or azide labeling reagent in water or DMSO
2. For each azide- or alkyne-modified cell or cell lysate sample, add the following reagents to a 1.5 mL microfuge tube, then vortex briefly to mix.
 - 50 µL cell or cell lysate sample
 - 50 µL PBS buffer
 - 50 µL of 5 mM corresponding dye azide (or dye alkyne) detection reagent in DMSO or water
3. Add 10 µL of 100 mM THPTA solution, vortex briefly to mix.
4. Add 10 µL of 20 mM CuSO₄ solution, vortex briefly to mix.
5. Add 10 µL of 300 mM sodium ascorbate solution to initiate the click reaction, vortex briefly to mix.
6. Protect the click reaction from light and allow it to incubate for 30 minutes at room temperature.
7. Cells or cell lysates are now click labeled and ready for downstream processing and/or analysis.

Appendix I. Chemical Properties of Tide Fluor™ Fluorescent Labeling Dyes

Tide Fluor™ dyes have improved labeling performance than the classic fluorescent labeling dyes such as FITC, TRITC, Texas Red®, Cy3, Cy5 and Cy7. They are the best affordable fluorescent dyes (alternative to Alexa Fluor® dyes) for labeling oligos and peptides without compromised performance. Each Tide Fluor™ dye is developed to match the spectral properties of a particular Alexa Fluor® or other labeling dyes (such as DyLight™ dyes).

| Labeling Dye | Cat# | Product Description | Reactivity | Adduct MW Calculation* |
|--------------|------|---|------------------------|------------------------|
| TF1 | 2236 | Tide Fluor™ 1 azide [TF1 azide] | Azide | + 301 |
| | 2237 | Tide Fluor™ 1 alkyne [TF1 alkyne] | Alkyne | + 270 |
| | 2238 | Tide Fluor™ 1 acid [TF1 acid] | NH ₂ and OH | + 215 |
| | 2239 | Tide Fluor™ 1 amine [TF1 amine] | CO ₂ H | + 257 |
| | 2242 | Tide Fluor™ 1 maleimide [TF1 maleimide] | SH | + 355 |
| | 2244 | Tide Fluor™ 1 succinimidyl ester [TF1 SE] | Aliphatic amine | + 215 |
| TF2 | 2245 | Tide Fluor™ 2 acid [TF2 acid] | NH ₂ and OH | + 469 |
| | 2246 | Tide Fluor™ 2 amine [TF2 amine] | CO ₂ H | + 398 |
| | 2247 | Tide Fluor™ 2 maleimide [TF2 maleimide] | SH | + 680 |
| | 2248 | Tide Fluor™ 2 succinimidyl ester [TF2 SE] | Aliphatic amine | + 469 |
| | 2252 | Tide Fluor™ 2 azide [TF2 azide] | Azide | + 555 |
| | 2253 | Tide Fluor™ 2 alkyne [TF2 alkyne] | Alkyne | + 524 |
| TF2WS | 2348 | Tide Fluor™ 2WS acid [TF2WS acid] | NH ₂ and OH | + 628 |
| | 2349 | Tide Fluor™ 2WS succinimidyl ester [TF2WS SE] | Aliphatic amine | + 628 |
| TF3 | 2254 | Tide Fluor™ 3 azide [TF3 azide] | Azide | + 526 |
| | 2255 | Tide Fluor™ 3 alkyne [TF3 alkyne] | Alkyne | + 495 |
| | 2268 | Tide Fluor™ 3 acid [TF3 acid] | NH ₂ and OH | + 440 |
| | 2269 | Tide Fluor™ 3 amine [TF3 amine] | CO ₂ H | + 496 |
| | 2270 | Tide Fluor™ 3 maleimide [TF3 maleimide] | SH | + 580 |
| | 2271 | Tide Fluor™ 3 succinimidyl ester [TF3 SE] | Aliphatic amine | + 440 |
| TF3WS | 2345 | Tide Fluor™ 3WS acid [TF3WS acid] | NH ₂ and OH | + 706 |
| | 2346 | Tide Fluor™ 3WS succinimidyl ester [TF3WS SE] | Aliphatic amine | + 706 |
| TF4 | 2285 | Tide Fluor™ 4 acid [TF4 acid] | NH ₂ and OH | + 544 |
| | 2286 | Tide Fluor™ 4 amine [TF4 amine] | CO ₂ H | + 586 |
| | 2287 | Tide Fluor™ 4 maleimide [TF4 maleimide] | SH | + 684 |
| | 2289 | Tide Fluor™ 4 succinimidyl ester [TF4 SE] | Aliphatic amine | + 544 |
| | 2300 | Tide Fluor™ 4 azide [TF4 azide] | Azide | + 630 |
| | 2301 | Tide Fluor™ 4 alkyne [TF4 alkyne] | Alkyne | + 599 |
| TF5WS | 2275 | Tide Fluor™ 5WS azide [TF5WS azide] | Azide | + 1078 |
| | 2276 | Tide Fluor™ 5WS alkyne [TF5WS alkyne] | Alkyne | + 787 |
| | 2278 | Tide Fluor™ 5WS acid [TF5WS acid] | NH ₂ and OH | + 732 |
| | 2279 | Tide Fluor™ 5WS amine [TF5WS amine] | CO ₂ H | + 774 |
| | 2280 | Tide Fluor™ 5WS maleimide [TF5WS maleimide] | SH | + 873 |
| | 2281 | Tide Fluor™ 5WS succinimidyl ester [TF5WS SE] | Aliphatic amine | + 732 |

[Continued on next page]

| | | | | |
|--------------|------|---|------------------------|--------|
| TF6WS | 2291 | Tide Fluor™ 6WS acid [TF6WS acid] | NH ₂ and OH | + 899 |
| | 2292 | Tide Fluor™ 6WS amine [TF6WS amine] | CO ₂ H | + 941 |
| | 2293 | Tide Fluor™ 6WS maleimide [TF6WS maleimide] | SH | + 1039 |
| | 2294 | Tide Fluor™ 6WS succinimidyl ester [TF6WS SE] | Aliphatic amine | + 899 |
| | 2302 | Tide Fluor™ 6WS azide [TF6WS azide] | Azide | + 1079 |
| | 2303 | Tide Fluor™ 6WS alkyne [TF6WS alkyne] | Alkyne | + 1048 |
| TF7WS | 2304 | Tide Fluor™ 7WS azide [TF7WS azide] | Azide | + 845 |
| | 2305 | Tide Fluor™ 7WS alkyne [TF7WS alkyne] | Alkyne | + 813 |
| | 2330 | Tide Fluor™ 7WS acid [TF7WS acid] | NH ₂ and OH | + 758 |
| | 2331 | Tide Fluor™ 7WS amine [TF7WS amine] | CO ₂ H | + 801 |
| | 2332 | Tide Fluor™ 7WS maleimide [TF7WS maleimide] | SH | + 899 |
| | 2333 | Tide Fluor™ 7WS succinimidyl ester [TF7WS SE] | Aliphatic amine | + 758 |
| TF8WS | 2306 | Tide Fluor™ 8WS azide [TF8WS azide] | Azide | + 1011 |
| | 2307 | Tide Fluor™ 8WS alkyne [TF8WS alkyne] | Alkyne | + 980 |
| | 2335 | Tide Fluor™ 8WS acid [TF8WS acid] | NH ₂ and OH | + 925 |
| | 2336 | Tide Fluor™ 8WS amine [TF8WS amine] | CO ₂ H | + 967 |
| | 2337 | Tide Fluor™ 8WS maleimide [TF8WS maleimide] | SH | + 1065 |
| | 2338 | Tide Fluor™ 8WS succinimidyl ester [TF8WS SE] | Aliphatic amine | + 925 |

* The molecular weight of the desired conjugate = the molecular weight of the free unlabeled molecule + the value listed in the table.

Appendix II. Spectral Properties of Tide Fluor™ Fluorescent Labeling Dyes

| Labeling Dye | Extinction Coefficient¹ (cm⁻¹M⁻¹) | Abs (nm) | Em (nm) | FQY² | CF at 260 nm³ | CF at 280 nm⁴ |
|---------------------|---|-----------------|----------------|------------------------|---------------------------------|---------------------------------|
| TF1 | 20,000 | 345 | 442 | 0.95 | 0.246 | 0.187 |
| TF2 | 75,000 | 500 | 527 | 0.90 | 0.288 | 0.201 |
| TF2WS | 75,000 | 502 | 525 | 0.90 | 0.211 | 0.091 |
| TF3 | 85,000 | 555 | 584 | 0.85 | 0.331 | 0.201 |
| TF3WS | 150,000 | 555 | 565 | 0.10 ⁵ | 0.079 | 0.079 |
| TF4 | 90,000 | 590 | 618 | 0.91 | 0.489 | 0.436 |
| TF5WS | 250,000 | 649 | 664 | 0.25 | 0.023 | 0.027 |
| TF6WS | 220,000 | 676 | 695 | 0.18 | 0.111 | 0.009 |
| TF7WS | 275,000 | 749 | 775 | 0.12 | 0.009 | 0.049 |
| TF8WS | 250,000 | 775 | 807 | 0.08 | 0.103 | 0.109 |

Notes: 1. Extinction Coefficient at their maximum absorption wavelength; 2. FQY = fluorescence quantum yield in aqueous buffer (pH 7.2); 3. CF at 260 nm is the correction factor used for eliminating the dye contribution to the absorbance at 260 nm (for oligo and nucleic acid labeling); 3. CF at 280 nm is the correction factor used for eliminating the dye contribution to the absorbance at 280 nm (for peptide and protein labeling); 5. Fluorescence intensity is significantly increased upon coupled to proteins or long peptides.

Appendix III. FRET Selection Guide of Tide Quencher™ Dyes

| Tide Fluor™ Donor | Ex(nm) | Em (nm) | Features and Benefits | Ordering Information |
|--|------------------|------------------|---|--|
| Tide Fluor™ 1 (TF1) | 345 nm | 442 nm | <i>Alternative to EDANS</i> <ul style="list-style-type: none"> • Much stronger absorption • Much stronger fluorescence • Less environment-sensitive | #2236 & 2237 (TF1 Click chemistry) #2238 (TF1 acid) #2239 (TF1 amine) #2242 (TF1 maleimide, SH-reactive) #2244 (TF1 SE, NH ₂ -reactive) |
| Tide Fluor™ 2 (TF2) Tide Fluor™ 2WS (TF2WS) | 500 nm 502 nm | 527 nm 525 nm | <i>Alternative to FAM, FITC and Alexa Fluor® 488</i> pH-insensitive fluorescence Photostable | #2245 (TF2 acid) & 2348 (TF2WS acid) #2246 (TF2 amine) #2247 (TF2 maleimide, SH-reactive) #2248 (TF2, SE) & #2249 (TF2WS SE) #2252 & 2253 (Click chemistry) |
| Tide Fluor™ 3 (TF3) Tide Fluor™ 3WS (TF3WS) | 555 nm 555 nm | 584 nm 565 nm | <i>Alternative to Cy3® and Alexa Fluor® 555</i> Strong fluorescence Photostable | #2254 & 2255 (TF3 Click chemistry) #2268 (TF3 acid) & 2345 (TF3WS acid) #2269 (TF3 amine) #2270 (TF3 maleimide, SH-reactive) #2271 (TF3 SE) & #2346 (TF3WS SE) |
| Tide Fluor™ 4 (TF4) | 590 nm | 618 nm | <i>Alternative to ROX, Texas Red® and Alexa Fluor® 594</i> Strong fluorescence Photostable | #2285 (TF4 acid) #2286 (TF4 amine) #2287 (TF4 maleimide, SH-reactive) #2289 (TF4 SE, NH ₂ -reactive) #2300 & 2301 (TF4 Click chemistry) |
| Tide Fluor™ 5WS (TF5WS) | 649 nm | 664 nm | <i>Alternative to Cy5® and Alexa Fluor® 647</i> Strong fluorescence Photostable | #2275 & 2276 (TF5WS Click chemistry) #2278 (TF5WS acid) #2279 (TF5WS amine) #2280 (TF5WS maleimide, SH-reactive) #2281 (TF5WS SE, NH ₂ -reactive) |
| Tide Fluor™ 6WS (TF6WS) | 676 nm | 695 nm | <i>Alternative to Cy5.5®, IRDye® 700 and Alexa Fluor® 680</i> Strong fluorescence Photostable | #2291 (TF6WS acid) #2292 (TF6WS amine) #2293 (TF6WS maleimide, SH-reactive) #2294 (TF6WS SE, NH ₂ -reactive) #2302 & 2303 (TF6WS Click chemistry) |
| Tide Fluor™ 7WS (TF7WS) | 749 nm | 775 nm | <i>Alternative to Cy7® and Alexa Fluor® 750</i> Strong fluorescence Photostable | #2304 & 2305 (TF7WS Click chemistry) #2330 (TF7WS acid) #2331 (TF7WS amine) #2332 (TF7WS maleimide, SH-reactive) #2333 (TF7WS SE, NH ₂ -reactive) |
| Tide Fluor™ 8WS (TF8WS) | 775 nm | 807 nm | <i>Alternative to IRDye® 800</i> Stronger fluorescence Higher Photostability | #2306 & 2307 (TF4 Click chemistry) #2335 (TF8WS acid) #2336 (TF8WS amine) #2337 (TF8WS maleimide, SH-reactive) #2338 (TF8WS SE, NH ₂ -reactive) |

*Texas Red® and Alexa Fluor® are the trademarks of Molecular Probes. CyDye, Cy3®, Cy5®, Cy5.5® and Cy7® are the trademarks of GE Health Care. IRDye® 700 and IRDye® 800 are the trademarks of Li-COR. Tide Fluor™ is the trademark of AAT Bioquest.

Appendix IV. HPLC Purification of Dye Oligonucleotide Conjugates

Ethanol Precipitation

Some commercial oligonucleotides often contain some interfering compounds, especially amines, such as triethylamine or Tris, and ammonium salts, we strongly recommend you to extract and precipitate the commercial oligo samples prior to initiating your labeling reaction. On the other hand, the labeling mixture contains labeled oligonucleotide, unlabeled oligonucleotide, hydrolyzed dye acid and unincorporated dye SE. The impurities of hydrolyzed dye acid and unincorporated dye SE resulted from the labeling reaction can be effectively removed by ethanol precipitation. The following protocol was optimized for the further purification of 0.1–1 mg commercial oligonucleotide sample that was purified by HPLC (3–30 A260 units).

- 1) Dissolve your target oligonucleotide in 100 μL of deionized water and extract three times with an equal volume of chloroform.
- 2) Precipitate the oligonucleotide by adding one-tenth volume (10 μL) of 3 M NaCl and two and a half volumes (250 μL) of cold absolute ethanol. Mix well and place at -20°C for 30 minutes.
- 3) Centrifuge the solution in a microcentrifuge at 10,000 to 15,000 g for 30 minutes.
- 4) Carefully remove the supernatant, rinse the pellet 1-3 times with cold 70% ethanol, and dry under a vacuum.
- 5) Dissolve the dry pellet in deionized water to achieve a final concentration of $>50 \mu\text{g}/\mu\text{L}$. This amine-modified oligonucleotide stock solution may be immediately used or stored frozen at $\leq -15^{\circ}\text{C}$.

Purification by HPLC

Labeled oligonucleotides can be purified by reverse-phase HPLC using a standard analytical C8 or C18 column using an analytical or semi-preparative HPLC instrument. The following protocol was optimized for the further purification of 0.1–1 mg labeled oligonucleotide (3–30 A260 units).

- 1) Dissolve the pellet from the ethanol precipitation in 0.1 M triethylammonium acetate (TEAA).
- 2) Load the dissolved pellet onto the column in 0.1 M TEAA and run a linear 5–95% acetonitrile gradient over 30 minutes.
Note 1: There will be peaks that correspond to the unlabeled oligonucleotide, the labeled oligonucleotide, and the free dye. The actual order and number of these peaks depends on the length of the oligonucleotide and the purity of the sample.
Note 2: To determine the identity of the peaks, monitor the absorbance at both 260 nm and at the absorbance maxima (λ_{max}) for the dye. For instruments with only one detector, two small samples should be run, each monitored at a different wavelength. Unlabeled oligonucleotide will show an absorbance at 260 nm only. Both the free dye and the labeled oligonucleotide will have absorbance at both 260 nm (A260 for oligo) and at the absorbance maximum of the dye (Amax for dye); The dye-labeled oligonucleotide will have a higher A260:Amax ratio than the dye or hydrolyzed dye.

Purification by Gel Electrophoresis

- 1) Pour a 0.5 mm-thick polyacrylamide slab gel.
Note: For oligonucleotides less than 25 bases in length, use 19% acrylamide, for oligonucleotides 25–40 bases, 15% acrylamide, and for oligonucleotides 40–100 bases, 12% acrylamide.
- 2) Resuspend the pellet from ethanol precipitation in 200 μL of 50% formamide, and incubate at 55°C for 5 minutes to disrupt any secondary structure.
- 3) Load the warmed oligonucleotide onto the gel and load an adjacent well with 50% formamide plus 0.05% bromophenol blue. The bromophenol blue will migrate at approximately the same rate as the oligonucleotide.
Note: You may need to use several wells.
- 4) Run the gel until the bromophenol blue indicator dye is two-thirds of the way down the gel.
- 5) Remove the gel from the glass plates and place on Saran Wrap.
- 6) Lay the gel on a fluorescent TLC plate.
- 7) Locate the labeled and unlabeled oligonucleotides by illumination with a handheld UV source.