

## Cell Meter™ Caspase 3/7 Activity Apoptosis Assay Kit \*Green Fluorescence\*

Catalog number: 22796  
Unit size: 200 Tests

Component	Storage	Amount
Component A: Caspase 3/7 Substrate (200X Stock Solution)	Freeze (<-15 °C), Minimize light exposure	2 vials (50 µL/vial)
Component B: Assay Buffer	Freeze (<-15 °C)	1 bottle (20 mL)

### OVERVIEW

Our Cell Meter™ assay kits are a set of tools for monitoring cell viability. There are a variety of parameters that can be used for monitoring cell viability. This particular kit is designed to monitor cell apoptosis through measuring Caspase 3 activation. Caspase 3 is widely accepted as a reliable indicator for cell apoptosis since the activation of caspase-3 (CPP32/apopain) is important for the initiation of apoptosis. Caspase 3 has substrate selectivity for the peptide sequence Asp-Glu-Val-Asp (DEVD). This kit uses Z-DEVD-Rh 110-DVED-Z as a fluorogenic indicator for caspase-3 activity. Cleavage of Rh 110 peptides by caspase 3 generates strongly fluorescent Rh 110 that is monitored fluorimetrically at 520-530 nm with excitation of 480-500 nm. The kit provides all the essential components with an optimized assay protocol. The assay is robust, and can be readily adapted for high-throughput assays. Using 100 µL of reagents per well in a 96-well format, this kit provides sufficient reagents to perform 200 assays. Using 25 µL of reagents per well in a 384-well format, this kit provides sufficient reagents to perform 800 assays.

### AT A GLANCE

#### Protocol summary

1. Prepare cells with test compounds (100 µL/well/96-well plate or 25 µL/well/384-well plate)
2. Add equal volume of Caspase 3/7 Substrate working solution
3. Incubate at room temperature for 1 hour
4. Monitor fluorescence intensity at Ex/Em = 490/525 nm (Cutoff = 515 nm)

**Important** Thaw one vial of each kit component at room temperature before starting the experiment.

### KEY PARAMETERS

Instrument:	Fluorescence microplate reader
Excitation:	490 nm
Emission:	525 nm
Cutoff:	515 nm
Recommended plate:	Black wall/clear bottom
Instrument specification(s):	Top/Bottom read mode

### PREPARATION OF WORKING SOLUTION

Add 50 µL of Caspase 3/7 Substrate (Component A) into 10 mL of Assay Buffer (Component B) and mix well to make Caspase 3/7 Substrate working solution.

**Note** Aliquot and store the unused Caspase 3/7 Substrate (Components A) and Assay Buffer (Component B) at -20 °C. Avoid repeated freeze/thaw cycles

### SAMPLE EXPERIMENTAL PROTOCOL

#### Cells preparation:

- **For adherent cells:** Plate cells overnight in growth medium at 20,000 cells/well/90 µL for a 96-well plate or 5,000cells/well/20mL for a 384-well plate.
- **For non-adherent cells:** Centrifuge the cells from the culture medium and then suspend the cell pellet in culture medium at 80,000 cells/well/90 µL for a

96-well poly-D lysine plate or 20,000 cells/well/20 µL for a 384-well poly-D lysine plate. Centrifuge the plate at 800 rpm for 2 minutes with brake off prior to the experiments.

**Note** Each cell line should be evaluated on an individual basis to determine the optimal cell density for apoptosis induction.


#### Sample Protocol:

1. Treat cells by adding 10 µL/well of 10X test compounds (96-well plate) or 5 µL/well of 5X test compounds (384-well plate) into PBS or the desired buffer. For blank wells (medium without the cells), add the same amount of compound buffer.
2. Incubate the cell plate in a 37°C, 5% CO<sub>2</sub>, incubator for a desired period of time (4 - 6 hours for Jurkat cells treated with camptothecin) to induce apoptosis.
3. Add 100 µL/well (96-well plate) or 25 µL/well (384-well plate) of Caspase 3/7 Substrate working solution.
4. Incubate the plate at room temperature for at least 1 hour, protected from light.

**Note** If desired, add 1 µL of the 1 mM Ac-DEVD-CHO caspase 3/7 inhibitor into selected samples 10 minutes before adding Caspase 3/7 working solution at room temperature to confirm the inhibition of the caspase 3/7-like activities.

5. Centrifuge cell plate (especially for the non-adherent cells) at 800 rpm for 2 minutes (brake off).
6. Monitor the fluorescence intensity with a fluorescence microplate reader at Ex/Em = 490/525 nm (Cutoff = 515 nm).

### EXAMPLE DATA ANALYSIS AND FIGURES



**Figure 1.** Detection of caspase 3/7 Activity in Jurkat cells. Jurkat cells were seeded on the same day at 80,000 cells/well/90 µL in a black wall/clear bottom 96-well costar plate. The cells were treated with or without 20 µM of camptothecin for 5 hours, and/or 5 µM caspase 3/7 inhibitor AC-DEVD-CHO for 10 minutes. The

caspace 3/7 working solution (100  $\mu$ L/well) was added and incubated at room temperature for 1 hour. The fluorescence intensity was measured at Ex/Em = 490/525 nm (Cutoff = 515 nm) using the NOVOstar instrument (BMG Labtech).

#### DISCLAIMER

AAT Bioquest provides high-quality reagents and materials for research use only. For proper handling of potentially hazardous chemicals, please consult the Safety Data Sheet (SDS) provided for the product. Chemical analysis and/or reverse engineering of any kit or its components is strictly prohibited without written permission from AAT Bioquest. Please call 408-733-1055 or email [info@aatbio.com](mailto:info@aatbio.com) if you have any questions.