

Cell Meter™ Autophagy Assay Kit *Green Fluorescence*

Catalog number: 23002
Unit size: 200 Tests

Component	Storage	Amount
Component A: 500X Autophagy Green™	Freeze (<-15 °C), Minimize light exposure	1 vial (50 µL)
Component B: Stain Buffer	Freeze (<-15 °C), Minimize light exposure	25 mL
Component C: Wash Buffer	Freeze (<-15 °C), Minimize light exposure	100 mL

OVERVIEW

This Cell Meter™ Autophagy Kit employs Autophagy Green™ as a specific autophagosome marker to analyze the activity of autophagy. The assay is optimized for direct detection of autophagy in both detached and attached cells. The kit provides all the essential components for the assay protocol. Cell Meter™ Autophagy Kit is suitable for fluorescence microscope, fluorescence microplate reader and flow cytometer. Autophagy Green™ has a large Stokes shift. Autophagy is an evolutionarily conserved degradation process that targets long-lived proteins, organelles, and other cytoplasmic components for degradation via the lysosomal pathway. The autophagy pathway is complementary to the action of the ubiquitin-proteasome pathway which typically degrades short-lived proteins. Activation of the autophagy pathway is required for multiple cellular roles, including survival during starvation, the clearance of intracellular components, development, and immunity. In the absence of stress, autophagy serves a house-keeping function, removing damaged organelles and cellular components preventing cytotoxic effects. Decreases and defects in autophagy have been implicated in multiple diseases, for example Huntingtons, Alzheimers, and Parkinsons. In terms of cancer development, autophagy seems to play multiple roles. Decreased or absent expression of certain autophagy proteins, such as Beclin-1 and Bif-1, increases tumor susceptibility in mice while the overexpression of these proteins can repress cancer cell growth. However, autophagy is critical for the survival of cancer cells within the nutrient poor and hypoxic core of solid tumors.

AT A GLANCE

Protocol summary

1. Prepare cells with your test compounds
2. Add Autophagy Green™ working solution
3. Incubate at 37°C for 15 - 60 minutes
4. Wash cells with Wash Buffer
5. Monitor the fluorescence increase at Ex/Em= 485/530 nm (Cutoff = 515 nm), fluorescence microscope with FITC filter set or flow cytometer with 530/30 nm filter (FITC channel)

Important Thaw all the components at room temperature before starting the experiment.

KEY PARAMETERS

Instrument:	Fluorescence microscope
Excitation:	FITC filter
Emission:	FITC filter
Recommended plate:	Black wall/clear bottom
Instrument:	Flow cytometer
Excitation:	488 nm laser
Emission:	530/30 nm filter
Instrument specification(s):	FITC channel
Instrument:	Fluorescence microplate reader
Excitation:	485 nm
Emission:	530 nm
Cutoff:	515 nm
Recommended plate:	Black wall/clear bottom
Instrument specification(s):	Bottom read mode

PREPARATION OF WORKING SOLUTION

Add 20 µL of Autophagy Green™ (Component A) to 10 mL of Stain Buffer (Component B) and mix well to make Autophagy Green™ working solution. Protect from light.

Note 20 µL of 500X Autophagy Green™ (Component A) is enough for one 96-well plate.

Note Aliquot and store unused 500X Autophagy Green™ at < -20°C. Protect from light and avoid repeated freeze-thaw cycles.

PREPARATION OF CELL SAMPLES

For guidelines on cell sample preparation, please visit <https://www.aatbio.com/resources/guides/cell-sample-preparation.html>

SAMPLE EXPERIMENTAL PROTOCOL

1. Culture cells to a density optimum for autophagy induction according to your specific induction protocol (about $1-2 \times 10^4$ cells/ well/96-well plate). At the same time, culture a non-induced negative control cell population at the same density as the induced population for every labeling condition.
2. Remove medium.
3. Add 100 µL/well (96-well plate) or 25 µL/well (384-well plate) of Autophagy Green™ working solution into each well.
4. Incubate the cells in a 37°C, 5% CO₂ incubator for 15 to 60 minutes.

Note The appropriate incubation time depends on the individual cell type and cell concentration used. Optimize the incubation time for each experiment.

5. Wash the cells with Wash Buffer (Component C) for 3 - 4 times, then add 100 µL Wash Buffer (Component C) to each well.

Note It is recommended to increase either the labeling concentration or the incubation time to allow the dye to accumulate if the cells do not appear to be sufficiently stained.

6. Monitor the fluorescence intensity with a fluorescence microplate reader at Ex/Em = 485/530 nm (Cutoff = 515 nm), a fluorescence microscope with FITC filter set or a flow cytometer with 530/30 nm filter (FITC channel).

EXAMPLE DATA ANALYSIS AND FIGURES

Control

Autophagy Treatment




Figure 1. Autophagy Green™ labeled vesicles were induced by starvation in HeLa cells. HeLa cells were incubated in a regular DMEM medium (Left: Control) or in 1X HBSS buffer with 5% serum (Right: Autophagy Treatment) for 16 hours. Both control and starved cells were incubated with Autophagy Green™ working solution for 20 minutes in a 37 °C, 5% CO₂ incubator, and then washed 3 times with wash buffer. Cells were imaged immediately under a fluorescence microscope with a FITC channel (green). Cell nuclei were stained with Hoechst 33342 (blue).

DISCLAIMER

AAT Bioquest provides high-quality reagents and materials for research use only. For proper handling of potentially hazardous chemicals, please consult the Safety Data Sheet (SDS) provided for the product. Chemical analysis and/or reverse engineering of any kit or its components is strictly prohibited without written permission from AAT Bioquest. Please call 408-733-1055 or email info@aatbio.com if you have any questions.