Screen QuestTM 10X Calcium Assay Buffer with Phenol Red PlusTM

Ordering Information:	Storage Conditions:	Instrument Platform:
Product Number: 36301 (10 plates)	Keep at -20 °C and protect from light	Florescence microplate readers

Introduction

Calcium flux assays are preferred methods in drug discovery for screening G protein coupled receptors (GPCR). Our Screen Quest™ 10X Calcium Assay Buffer with Phenol Red Plus™ contains our water soluble and heat stable probenecid which inhibits the activities of drug-efflux pumps. It can be used to prevent florescent dyes (such as Indo-1 AM, Fura-2 AM, Fluo-3 AM, Fluo-4 AM, Fluo-8 AM, Rhod-2 AM and Rhod-4 AM) from leaking out of cells.

Kit Component

Component	Amount
Screen Quest [™] 10X Calcium Assay Buffer with Phenol Red Plus [™]	1 bottle (10 mL)

Protocol (for one plate)

- 1. Thaw the bottle at room temperature before use.
- 2. <u>Make 1X Screen Quest TM calcium assay buffer:</u> Add 1 mL of 10X Screen Quest TM calcium assay buffer with Phenol Red PlusTM (Cat. # 36301) to 9 mL of HHBS (1X Hank's with 20 mM Hepes buffer, pH 7.0), and mix them well.
 - Note: 10 mL of 1X assay buffer is enough for one plate. Aliquot and store unused 10X assay buffer at \leq -20 °C. Protect from light and avoid repeated freeze-thaw cycles.
- 3. Make 2X dye-loading solution for one cell plate: Add DMSO reconstituted fluorescent calcium dyes (such as Indo-1 AM, Fura-2 AM, Fluo-3 AM, Fluo-4 AM and Fluo-8 AM, Rhod-2 AM and Rhod-4 AM) into 10 mL of 1X Screen Quest TM calcium assay buffer (from Step 2), to make the final well dye concentration 2X of the desired concentration, and mix them well. The working solution is stable for at least 2 hours at room temperature.
- 4. To the microplate well add 2X dye-loading solution (from Step 3) which is the same volume as the cell culture medium (e.g., 100 uL/well/96-well or 25 μL/well/384-well).
- 5. Incubate the cells in a 37 °C, 5% CO₂ incubator for about 1 hours, or as desired.
- 6. Prepare the compound plate with HHBS or your desired buffer.
- 7. Run the calcium flux assay.

.