

## Screen Quest™ Fluorimetric Glucose Uptake Assay Kit

Catalog number: 36500, 36501  
Unit size: 100 Tests, 500 Tests

Component	Storage	Amount	
		Cat No. 36500	Cat No. 36501
Component A: 2-Deoxyglucose (2-DG, 10mM )	Freeze (<-15 °C), Minimize light exposure	1 mL	5 mL
Component B: Glucose Uptake Buffer	Freeze (<-15 °C), Minimize light exposure	10 mL	50 mL
Component C: Acidic Lysis Buffer	Freeze (<-15 °C), Minimize light exposure	2.5 mL	12.5 mL
Component D: Neutralization Buffer	Freeze (<-15 °C), Minimize light exposure	2.5 mL	12.5 mL
Component E: Enzyme Probe	Freeze (<-15 °C), Minimize light exposure	1 bottle (lyophilized powder)	5 bottles (lyophilized powder)
Component F: Assay Buffer	Freeze (<-15 °C)	1 bottle (5 mL)	1 bottle (25 mL)
Component G: NADP	Freeze (<-15 °C), Minimize light exposure	1 vial	5 vials
Component H: 5x KRPH Buffer	Freeze (<-15 °C), Minimize light exposure	20 mL	100 mL

### OVERVIEW

Glucose transport systems are responsible for transporting glucose across cell membranes. Measuring uptake of 2-deoxyglucose (2-DG), a glucose analog, in tissues and cells is widely accepted as a reliable method to estimate the amount of glucose uptake and to investigate the regulation of glucose metabolism and mechanism of insulin resistance. The 2-DG uptake is commonly determined by using non-metabolized 2-DG labeled with tritium or C14. However, routine use of a radiolabelled probe is costly and requires a tedious special handling procedure. AAT Bioquest's Screen Quest™ Fluorimetric Glucose Uptake Assay Kit provides a sensitive and non-radioactive assay in cells. In this assay 2-DG is taken up by glucose transporters, and metabolized to 2-DG-6-phosphate (2-DG6P). The non-metabolizable 2-DG6P accumulates in the cells, and is proportional to glucose uptake by cells. The accumulated 2-DG6P is enzymatically oxidized and generates NADPH, which is specifically monitored by a fluorogenic NADPH sensor. The signal can be read by a fluorescence microplate reader.

### AT A GLANCE

#### Protocol summary

1. Plate cells and treat the cells as desired
2. Add 10 µL/well 2-DG and incubate at 37°C for 20-40 minutes
3. Wash cells and lyse cells
4. Add 50 µL/well of 2-DG Uptake Assay working solution
5. Incubate at RT for 30 to 120 minutes
6. Monitor fluorescence increase at Ex/Em = 540/590 nm

**Important** Thaw all the kit components at room temperature before starting the experiment. The above mentioned protocol is for one 96-well plate.

### KEY PARAMETERS

Instrument:	Fluorescence microplate reader
Excitation:	540 nm
Emission:	590 nm
Cutoff:	570 nm
Recommended plate:	Solid black

### PREPARATION OF STOCK SOLUTIONS

*Unless otherwise noted, all unused stock solutions should be divided into single-use aliquots and stored at -20 °C after preparation. Avoid repeated freeze-thaw cycles.*

#### KRPH Buffer stock solution (1X):

Add 20 mL of KRPH Buffer (5X) (Component H) to 80 mL of deionized water and

mix them well.

**Note** 50 mL volume of 1× KRPH Buffer is enough for approximately one 96-well plate. Prepare the needed volume proportionally. Store the unused 1×KRPH at 4°C or -20°C.

### PREPARATION OF WORKING SOLUTION

#### 1. NADP working solution:

Add 100 µL of H<sub>2</sub>O into the vial of NADP (Component G) and mix them well.

#### 2. Enzyme Probe working solution:

Add 5 mL of Assay Buffer (Component F) into the bottle of Enzyme Probe (Component E) and mix them well.

#### 3. 2-DG Uptake Assay working solution:

Add 100 µL of NADP working solution into the Enzyme Probe working solution and mix them well.

**Note** These quantities are good for one 96-well plate.

### SAMPLE EXPERIMENTAL PROTOCOL

**Important:** This protocol can be used as guidelines to culture 3T3-L1 adipocytes for 2-DG uptake.


1. Plate cells in growth medium at 50,000-80,000 cells/well/100 µL/96-well or 12,500-20,000 cells/well/25 µL/384-well black wall/clear bottom cell culture Poly-D lysine plate for 4-6 hours before experiment.
2. Remove the cell plate from the incubator, aspirate the medium from the wells, and deprive the cells with 100 µL/well (96 well-plate) or 25 µL/well (384 well-plate) serum free medium. Incubate the cells at 37°C, 5% CO<sub>2</sub> incubator for 6 hours to overnight.
3. Remove the cell plate from the incubator, aspirate the medium from the wells, and gently wash the cells twice with 100 µL/well 1× KRPH buffer.
4. Add 90 µL/well Glucose Uptake Buffer (Component B) and incubate the cells at 37°C, 5% CO<sub>2</sub> incubator for 1 hour.
5. Stimulate with or without insulin or compound of test for 20 mins. Add 10 µL/well of the 10× insulin solution to a final concentration of 1 µM or 10× compound solution of test. And also add 10 µL insulin vehicle buffer or compound vehicle buffer to the untreated wells as control, and incubate at 37°C, 5% CO<sub>2</sub> incubator for 20 mins.

- For glucose uptake inhibition study, add 10× Phloretin to a final concentration of 200 μM or inhibitors of test, and incubate at 37°C, 5% CO<sub>2</sub> for 2-5 min.

**Note** 10 mL inhibitor vehicle buffer is suggested to be added to both the insulin treated and untreated wells as control. Phloretin treated cells can be used as positive control.

- Add 10 μL/well 2-DG solution (Component A) to each well, and incubate at 37°C, 5% CO<sub>2</sub> incubator for 20-40 min. For negative controls, leave some wells untreated with insulin, inhibitor and 2-DG.
- After treatment, remove solution in each well and gently wash cells 3 times, 100 μL/well with KRPH Buffer to remove the extra 2-DG from the solution. Remove KRPH Buffer from the wells.
- Add 25 μL/well Acidic Lysis Buffer (Component C) to each well and incubate at 37°C for 20 mins to lyse the cells. And the 2DG uptake assay mixture could be prepared in the meantime.
- Add 25 μL/well Neutralization Buffer (Component D) to each well, mix thoroughly, leave at room temperature for 5-10 minutes to neutralize the cell lysate.
- Add 50 μL of 2DG Uptake Assay working solution to each well of 2DG6P standard (Not provided) or cell lysate.
- Incubate the reaction at room temperature for 30 minutes to 2 hours, protected from light.
- Monitor the fluorescence increase with a fluorescence plate reader at Ex/Em = 530-570/590-600 nm (optimal Ex/Em = 540/590 nm, cut off at 570 nm).

#### EXAMPLE DATA ANALYSIS AND FIGURES



**Figure 1.** Measurement of 2DG uptake in differentiated 3T3-L1 adipocytes and 3T3-L1 fibroblasts. Assays were performed with Screen Quest™ Fluorimetric Glucose Uptake Assay Kit in a black wall/clear bottom cell culture Poly-D lysine plate using a Gemini (Molecular Devices) microplate reader. (A: Negative Control, no insulin no 2-DG treatment. B: 2DG uptake in the absence of insulin. C: 2DG uptake in the presence of 1mM insulin. D: 2DG uptake in the presence of 1mM insulin and 200 mM phloretin. E: 2DG uptake in the presence of insulin 1mM and 5mM D-Glucose.)

#### DISCLAIMER

AAT Bioquest provides high-quality reagents and materials for research use only. For proper handling of potentially hazardous chemicals, please consult the Safety Data Sheet (SDS) provided for the product. Chemical analysis and/or reverse engineering of any kit or its components is strictly prohibited without written permission from AAT Bioquest. Please call 408-733-1055 or email [info@aatbio.com](mailto:info@aatbio.com) if you have any questions.