

iFluor™ 860 acid

Catalog number: 1405

Unit size: 1 mg

Product Details

Storage Conditions Freeze (<-15 °C), Minimize light exposure

Expiration Date 12 months upon receiving

Chemical Properties

Appearance Dark green solid

Molecular Weight 1403.41

Soluble In DMSO

Spectral Properties

Excitation Wavelength 853 nm

Emission Wavelength 878 nm

Applications

In vivo fluorescence imaging uses a sensitive camera to detect fluorescence emission from fluorophores in whole-body living small animals. To overcome the photon attenuation in living tissue, fluorophores with long emission at the infrared (IR) region are generally preferred. Recent advances in imaging strategies and reporter techniques for in vivo fluorescence imaging include novel approaches to improve the specificity and affinity of the probes and to modulate and amplify the signal at target sites for enhanced sensitivity. Further emerging developments are aiming to achieve high-resolution, multimodality and lifetime-based in vivo fluorescence imaging. Our iFluor™ 860 is designed to label proteins and other biomolecules with infrared fluorescence. Conjugates prepared with iFluor™ 860 have the excitation and emission in the IR range. iFluor™ 860 dye emission is well separated from commonly used far-red fluorophores such as Cy5, Cy7 or allophycocyanin (APC), facilitating multicolor analysis. This fluorophore is also useful for small animal in-vivo imaging applications or for other imaging applications that require IR detections.