

ReadiLink™ Rapid XFD488 Antibody Labeling Kit *XFD488 Same Structure to Alexa Fluor™ 488*

Catalog number: 1275
Unit size: 2x50 ug Labelings

Product Details

Expiration Date 6 months upon receiving

Spectral Properties

Excitation Wavelength 499 nm

Emission Wavelength 520 nm

Applications

XFD488 is manufactured by AAT Bioquest, and it has the same chemical structure of Alexa Fluor® 488 (Alexa Fluor® is the trademark of ThermoFisher). Readilink™ protein labeling technology is the most robust and convenient tool for preparing fluorescent antibody conjugates used for fluorescence imaging and flow cytometry applications. ReadiLink™ Rapid XFD488 Antibody Labeling Kit provides the most convenient tool for making XFD488-labeled antibody conjugates. XFD488 dye used in this ReadiLink™ kit is reasonably stable and shows good reactivity and selectivity with protein amino groups. The kit has all the essential components for labeling ~2x50 ug antibody. Each of the two vials of XFD488 dye provided in the kit is optimized for labeling ~50 µg antibody. ReadiLink™ Rapid XFD488 Antibody Labeling Kit requires minimal hands-on time. The prepared XFD488 conjugates (with the kit) are ready to use for fluorescence imaging and flow cytometry applications without further purifications needed. It provides the most convenient method to prepare XFD488-labeled antibodies.