

**XFD350 goat anti-rabbit IgG (H+L) *Cross
Adsorbed, XFD350 Same Structure to Alexa
Fluor™ 350***Catalog number: 16395
Unit size: 1 mg**Product Details**

Storage Conditions	2-6°C and kept from light. To extend the shelf-life of this product, add an equal volume of glycerol to make a final concentration of approximately 50% glycerol and store at -20°C.
Expiration Date	12 months upon receiving
Concentration	1 mg/mL
Formulation	PBS, 2 mg/mL BSA

Unit Details

Unit	16395 (1 mg)
Reconstitution Volume	1 mL ddH ₂ O

Antibody Properties

Species Reactivity	Rabbit
Class	Secondary
Clonality	Polyclonal
Host	Goat

Chemical Properties

Molecular Weight	~150000
------------------	---------

Biological Properties

Stabilizer	None
Appearance	Off-white solid
Preparation	Goat anti-rabbit IgG (H+L) is produced in goat with pooled total rabbit IgG, and affinity purified with rabbit IgG coupled beads. The purified IgG has a minimal cross-reaction to human, horse, mouse, human and bovine IgG. The antibody is conjugated with XFD350 under optimal condition.
Application	Flow Cytometry (FACS), ELISA, HC, Western Blot
Soluble In	Water

Spectral Properties

Conjugate	XFD350
Excitation Wavelength	343 nm
Emission Wavelength	441 nm

Applications

XFD350 is manufactured by AAT Bioquest, and it has the same chemical structure of Alexa Fluor® 350 (Alexa Fluor® is the trademark of ThermoFisher). Our goat anti-rabbit IgG whole antibodies have been cross-adsorbed against human IgG and human serum prior to conjugation to minimize cross-reactivity. This XFD350 labeled-goat anti-rabbit IgG conjugate is prepared by the reaction of cross-adsorbed goat anti-rabbit IgG whole antibody with XFD350 NHS ester, which has the same molecule to Alexa Fluor® 350 NHS ester. Each conjugate has typically 4-6 fluorophores per IgG molecule. Fluorescent secondary antibody conjugates are useful in the detection, sorting, or purification of its specified target and ideal for fluorescence microscopy and confocal laser scanning microscopy, flow cytometry, and fluorescent western detection.