

iFluor™ 568 goat anti-mouse IgG (H+L)

Catalog number: 16462, 16470

Unit size: 200 ug, 1 mg

Product Details

Storage Conditions 2-6°C and kept from light. To extend the shelf-life of this product, add an equal volume of

glycerol to make a final concentration of approximately 50% glycerol and store at -20°C.

Expiration Date 12 months upon receiving

Concentration 1 mg/mL

Formulation PBS, 2 mg/mL BSA

Unit Details

Unit 16462 (200 ug) 16470 (1 mg)

Reconstitution Volume 200 uL ddH₂O 1 mL ddH₂O

Antibody Properties

Species Reactivity Mouse

Class Secondary

Clonality Polyclonal

Host Goat

Chemical Properties

Molecular Weight ~150000

Biological Properties

Stabilizer None

Appearance Purple solid

Preparation Goat anti-mouse IgG (H+L) is produced in goat with pooled total mouse IgG, and affinity purified

with mouse IgG coupled beads. The antibody is conjugated with iFluor™ 568 under optimal

condition.

Application Immunofluorescence (IF), Flow Cytometry (FACS)

Soluble In Water

Spectral Properties

Conjugate iFluor™ 568

Excitation Wavelength 568 nm

Applications

iFluor™ 568 is a bright red fluorescent dye. iFluor™ 568-labeled anti-IgG conjugates exhibit bright fluorescence signal and good photostability. Used for stable signal generation in imaging and flow cytometry, the fluorescence intensity of iFluor™ 568 conjugates is pH-insensitive from pH 4 to pH 11. The iFluor™ 568-labeled antibody conjugates can be well excited with Krypton ion laser (~568 nm). iFluor™ 568 family has the spectral properties essentially identical to those of Alexa Fluor® 568. Under the same conditions we tested, iFluor™ 568 antibody conjugates are brighter and more photostable than the corresponding Alexa Fluor® 568. These spectral and labeling characteristics make the iFluor™ 568 dye family a superior alternative to Alexa Fluor® 568. In addition, iFluor™ 568 secondary antibody conjugates give higher signal/background ratios than the corresponding Alexa Fluor® 568-labeled conjugates.