

**iFluor™ 546 goat anti-rabbit IgG (H+L)**Catalog number: 16618  
Unit size: 200 ug**Product Details**

---

Storage Conditions	2-6°C and kept from light. To extend the shelf-life of this product, add an equal volume of glycerol to make a final concentration of approximately 50% glycerol and store at -20°C.
Expiration Date	12 months upon receiving
Concentration	1 mg/mL
Formulation	PBS, 2 mg/mL BSA

**Unit Details**

---

Unit	16618 (200 ug)
Reconstitution Volume	200 uL ddH <sub>2</sub> O

**Antibody Properties**

---

Species Reactivity	Rabbit
Class	Secondary
Clonality	Polyclonal
Host	Goat

**Chemical Properties**

---

Molecular Weight	~150000
------------------	---------

**Biological Properties**

---

Stabilizer	None
Appearance	Red solid
Preparation	Goat anti-rabbit IgG (H+L) is produced in goat with pooled total rabbit IgG, and affinity purified with rabbit IgG coupled beads. The antibody is conjugated with iFluor™ 546 under optimal condition.
Application	Immunofluorescence (IF), Flow Cytometry (FACS)
Soluble In	Water

**Spectral Properties**

---

Conjugate	iFluor™ 546
-----------	-------------

Excitation Wavelength	541 nm
Emission Wavelength	557 nm

## Applications

---

iFluor™ 546 is a bright orange fluorescent dye. iFluor™ 546-labeled anti-IgG conjugates exhibit bright fluorescence signal and good photostability. Used for stable signal generation in imaging and flow cytometry, the fluorescence intensity of iFluor™ 546 conjugates is pH-insensitive from pH 4 to pH 11. The iFluor™ 546-labeled antibody conjugates can be well excited with either Nd:YAG laser (~532 nm) or Helium-Neon laser (~543 nm). iFluor™ 546 family has the spectral properties essentially identical to those of Alexa Fluor® 546. Under the same conditions we tested, iFluor™ 546 antibody conjugates are brighter and more photostable than the corresponding Alexa Fluor® 546. These spectral and labeling characteristics make the iFluor™ 546 dye family a superior alternative to Alexa Fluor® 546. In addition, iFluor™ 546 secondary antibody conjugates give higher signal/background ratios than the corresponding Alexa Fluor® 546-labeled conjugates.