

iFluor™ 568 goat anti-rabbit IgG (H+L)Catalog number: 16622, 16804
Unit size: 200 ug, 1 mg**Product Details**

Storage Conditions	2-6°C and kept from light. To extend the shelf-life of this product, add an equal volume of glycerol to make a final concentration of approximately 50% glycerol and store at -20°C.
Expiration Date	12 months upon receiving
Concentration	1 mg/mL
Formulation	PBS, 2 mg/mL BSA

Unit Details

Unit	16622 (200 ug)	16804 (1 mg)
Reconstitution Volume	200 uL ddH ₂ O	1 mL ddH ₂ O

Antibody Properties

Species Reactivity	Rabbit
Class	Secondary
Clonality	Polyclonal
Host	Goat

Chemical Properties

Molecular Weight	~150000
------------------	---------

Biological Properties

Stabilizer	None
Appearance	Purple solid
Preparation	Goat anti-rabbit IgG (H+L) is produced in goat with pooled total rabbit IgG, and affinity purified with rabbit IgG coupled beads. The antibody is conjugated with iFluor™ 568 under optimal condition.
Application	Immunofluorescence (IF), Flow Cytometry (FACS)
Soluble In	Water

Spectral Properties

Conjugate	iFluor™ 568
Excitation Wavelength	568 nm

Applications

iFluor™ 568 is a bright red fluorescent dye. iFluor™ 568-labeled anti-IgG conjugates exhibit bright fluorescence signal and good photostability. Used for stable signal generation in imaging and flow cytometry, the fluorescence intensity of iFluor™ 568 conjugates is pH-insensitive from pH 4 to pH 11. The iFluor™ 568-labeled antibody conjugates can be well excited with Krypton ion laser (~568 nm). iFluor™ 568 family has the spectral properties essentially identical to those of Alexa Fluor® 568. Under the same conditions we tested, iFluor™ 568 antibody conjugates are brighter and more photostable than the corresponding Alexa Fluor® 568. These spectral and labeling characteristics make the iFluor™ 568 dye family a superior alternative to Alexa Fluor® 568. In addition, iFluor™ 568 secondary antibody conjugates give higher signal/background ratios than the corresponding Alexa Fluor® 568-labeled conjugates.