

iFluor™ 594 goat anti-rabbit IgG (H+L)

Catalog number: 16628, 16806

Unit size: 200 ug, 1 mg

Product Details

Storage Conditions	2-6°C and kept from light. To extend the shelf-life of this product, add an equal volume of glycerol to make a final concentration of approximately 50% glycerol and store at -20°C.
Expiration Date	12 months upon receiving
Concentration	1 mg/mL
Formulation	PBS, 2 mg/mL BSA

Unit Details

Unit	16628 (200 ug)	16806 (1 mg)
Reconstitution Volume	200 uL ddH ₂ O	1 mL ddH ₂ O

Antibody Properties

Species Reactivity	Rabbit
Class	Secondary
Clonality	Polyclonal
Host	Goat

Chemical Properties

Molecular Weight	~150000
------------------	---------

Biological Properties

Stabilizer	None
Appearance	Purple solid
Preparation	Goat anti-rabbit IgG (H+L) is produced in goat with pooled total rabbit IgG, and affinity purified with rabbit IgG coupled beads. The antibody is conjugated with iFluor™ 594 under optimal condition.
Application	Flow Cytometry (FACS), ELISA, HC, Western Blot
Soluble In	Water

Spectral Properties

Conjugate	iFluor™ 594
Excitation Wavelength	588 nm

Applications

AAT Bioquest's iFluor™ dyes are optimized for labeling proteins, in particular, antibodies. These dyes are bright, photostable and have minimal quenching on proteins. They can be well excited by the major laser lines of fluorescence instruments (e.g., 350, 405, 488, 555 and 633 nm). iFluor™ 594 goat anti-rabbit IgG (H+L) conjugate has fluorescence excitation and emission maxima of ~592 nm and ~614 nm respectively. These spectral characteristics make them an excellent alternative to Alexa Fluor® 594 goat anti-rabbit IgG (H+L) conjugate (Alexa Fluor® is the trademark of Invitrogen).