

MagaDye™ 535-ddGTP

Catalog number: 17063

Unit size: 5 nmoles

Product Details

Storage Conditions Freeze (<-15 °C), Minimize light exposure

Expiration Date 12 months upon receiving

Chemical Properties

Appearance Liquid

Soluble In Water

Applications

Sanger sequencing, also known as the chain termination method, is a technique for DNA sequencing based upon the selective incorporation of chain-terminating dideoxynucleotides (ddNTPs) by DNA polymerase. It was developed by Frederick Sanger and colleagues in 1977. Although the newer NGS technologies are becoming common in clinical research labs due to their higher throughput capabilities and lower costs per sample, Sanger sequencing with 99.99% accuracy is still the “gold standard” for clinical research sequencing. The four distinct fluorescent ddNTPs (labeled with BigDye®, BigDye® is the trademark of ThermoFisher) are the critical components for performing Sanger sequencing. MagaDye™ 535-ddGTP is equivalent to BigDye dR6G color with almost identical spectra.