

Product Information Sheet

Ordering Information

50006
3-Amino-2-oxazolidine (AOZ)-HRP Conjugate
mg
200uL in 100mM Sodium Phosphate Buffer, pH 6.0
2mg/ml BSA, 3% trehalose, and 0.01% thimerosal after reconstitution
Liquid
Store at Refrigerator (2-8 °C)
20 21

Application Notes

AOZ-HRP conjugates have been characterized with our ABTS (#11001), TMB (#11003) and ADHP (#11000) peroxidase substrates. A typical of 1:5,000 to 20,000 dilution were used for ELISA application

*The suggested working dilution is provided as a guide only. It is recommended that the users titrates the product for use in their tests using proper positive and negative controls.